
Latex-Free CPR Bag
Disposable Manual Resuscitator
DIRECTIONS FOR USE
Adult - Small Adult - Child - Infant
Non-Sterile • Single-Patient-Use
Notice: Should be used only by personnel trained in cardiopulmonary resuscitation.


- PRECAUTIONS:**
- Resuscitator should be used only by personnel trained in cardiopulmonary resuscitation.
 - Never wait to begin mouth-to-mask resuscitation if a manual resuscitator is not immediately available or cannot be used effectively. (See your department manual for accepted patient resuscitation procedure.)
 - If used with supplemental oxygen, do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable chemicals.
 - Should not be used in toxic or hazardous atmospheres.
 - Infant and Child units are equipped with a pressure-limiting device which opens at approximately 40 cm H₂O. However, an abrupt, high-volume inspiratory effort may cause the unit to exceed this level.
 - Verify proper function of the resuscitator and clear patient airway by monitoring for the following:
 - That the patient is being ventilated, as indicated by rise and fall of chest. The use of an airway pressure manometer is recommended.
 - Proper CPR Bag valve action. (Refer to #3. Prior to using, TEST THE RESUSCITATOR.)
 - That the mask and valve are free from obstruction. To clear valve obstructions, squeeze and shake the bag briskly or rinse with water.
 - Do not sterilize the CPR Bag. This product is designed for single-patient-use only and parts are not designed for re-use or to be interchanged with the permanent LAERDAL SILICONE RESUSCITATOR.

- DIRECTIONS FOR USE:**
- Remove the resuscitator from the outer protective poly bag. Expand the CPR Bag to its operating position.
 - Inspect the unit to be sure the system is complete.
 - Prior to using, TEST THE RESUSCITATOR:

a. Compress the ventilation bag with one hand, then occlude patient valve outlet with your other hand. Release the grip on the bag. Rapid bag re-expansion confirms efficient air intake.


- b. Remove patient valve, close the neck opening and try to compress the bag. If the bag cannot be compressed with reasonable force, or if bag compression forces the air out between your hand and neck of the bag, the intake valve is efficiently preventing backward escape of air.
- c. Attach the patient valve to the bag. Place a test lung over the valve connector (patient side). Test the resuscitator by squeezing and releasing the resuscitator bag with one hand as rapidly as possible for at least 10 breaths. Fully compress the resuscitator bag with each squeeze. This should fill the test lung and confirm that the patient valve efficiently directs inspiration air to the patient. If the test lung continues to inflate because of the stacking of breaths, suspect a faulty exhalation valve, and remove the resuscitator from service.
- d. Compress the filled test lung. Air should vent to the atmosphere through the deflector and not return to the ventilation bag. When using unit with a pressure-limiting device, test for proper function by occluding patient valve outlet and compressing the bag to verify opening of the pressure-limiting device.

Important: If the resuscitator does not pass all criteria as described under step #3 "Test the Resuscitator" take immediate action by replacing resuscitator.

4. **INSPIRATION PRESSURE LIMITING DEVICE:** The infant and child resuscitators may feature a patient valve with a special pressure-limiting device mounted on the upper valve housing. If inspiration meets with pulmonary resistance, venting will occur, limiting pressure to 40cmH₂O, thereby reducing the risk of stomach distension. A hissing sound can be heard as the device opens. **THE CHILD AND INFANT VALVE PROVIDES A BUILT-IN PORT FOR PRESSURE MONITORING. RECOMMENDED WHEN RESUSCITATING INFANTS AND SMALL CHILDREN, REMOVE THE CAP AND ATTACH YOUR MONITORING DEVICE. REATCHAT CAP WHEN NOT USING.**

- Note:** When higher inspiration pressures are necessary, the pressure limiting device can be closed with the tip of the index finger or by engaging the lock-off clip while squeezing the bag. The use of an airway pressure manometer is recommended when doing this procedure.
- If resuscitating where high concentrations of oxygen are needed, attach oxygen tubing adapter to proper oxygen source, such as a flow meter or adjustable oxygen regulator.
 - When a bag reservoir is being used, adjust oxygen flow as needed so the bag will remain fully or partially inflated during use, whichever is preferred. When an aerosol reservoir is being used, a 15 liter flow is recommended.
 - When using a mask attached to the CPR Bag, be sure you have a tight, secure fit.
 - When using an endotracheal tube or tracheotomy tube adapter, remove mask and attach the patient valve outlet directly to the tube. The swivel connection provides a 15mm I.D. port for this purpose.
 - The adjusting Valve features a built-in PEEP Adapter Port. To attach the PEEP Valve to the Port, remove the exhalation deflector and add the PEEP Valve with a slight twist as you push together the parts to assure a tight, secure fit. When adjusting the MERCURY MEDICAL® PEEP Valve, connect a manometer in line with the patient breathing system and read the manometer during adjustment. Rotate adjustment cap clockwise to increase PEEP or counterclockwise to decrease PEEP. The adjustment range is 0 - 20 cm H₂O (mbar).

Warning: The MERCURY MEDICAL® PEEP Valve should be used only by experienced persons who are aware that Positive End Expiratory Pressure (PEEP) may have an adverse effect on the patient's circulation. When using PEEP, the patient's condition should be monitored and PEEP pressure verified with a manometer.

- OPERATING INSTRUCTIONS:**
- Clear mouth of foreign matter.
 - Position patient for open airway.
 - Apply mask firmly to face to achieve a tight seal. If patient is intubated, attach the patient valve connector to the tube adapter.
 - Squeeze and release bag quickly, allowing enough time between squeezes for the patient to exhale and the bag to re-expand.
 - Observe rise and fall of patient's chest and listen for air flow from patient's valve as patient exhales.
 - Verwijder vreemde deeltjes uit de mond.
 - Plaats de patiënt zodanig dat de luchtweg open blijft.
 - Zet het masker stevig op het gezicht en voer een goede afsluiting toe. Zorg dat de patiënt voldoende tijd heeft om uit te ademen, en dat de patiëntklep weer sluit.
 - Kijk naar het omhoog en omlaag gaan van de borstkas van de patiënt en luister naar de lucht die eruit of juist in de patiëntklep terechtkomt.

Important: If these do not occur, patient's airway or patient valve may be blocked. Take immediate action by replacing the resuscitator or use an alternative procedure appropriate for situation.

SPECIFICATIONS AND PERFORMANCE:
The MERCURY MEDICAL® CPR Bag was designed to meet ASTM 920-93 and ISO 8382 requirements.
Note: Performance values given are achievable under test conditions but may vary during actual use. Information on test methods is available from MERCURY MEDICAL®. Performance specifications verified by independent test laboratory.

| OXIGEN CONCENTRATION: | | LPM | | Rate | | Tidal Volume | | w/reservoir | | w/o reservoir | |
|--------------------------|-----|--------|--------------|-------------|-----|--------------|--------|-------------|-----|---------------|-----|
| Adult & Small Adult Unit | 5 | 12 | 500 ml | 68% | 34% | 10 | 12 | 500 ml | 92% | 47% | 70% |
| | 10 | 12 | 500 ml | 100% | 70% | 12 | 12 | 500 ml | 95% | 96% | 96% |
| | 15 | 12 | 500 ml | 97% | | 15 | 12 | 500 ml | 97% | | |
| Child Unit | LPM | Rate | Tidal Volume | w/reservoir | 10 | 10 | 12 | 500 ml | 92% | | |
| | 20 | 250 ml | 250 ml | 85% | 12 | 12 | 500 ml | 95% | | | |
| | 14 | 12 | 500 ml | 96% | 14 | 12 | 500 ml | 96% | | | |
| | 15 | 12 | 500 ml | 97% | 15 | 12 | 500 ml | 97% | | | |
| Infant Unit | LPM | Rate | Tidal Volume | w/reservoir | 10 | 10 | 12 | 500 ml | 92% | | |
| | 4 | 30 | 40 ml | 100% | 85% | 12 | 12 | 500 ml | 95% | | |

| MAXIMUM STROKE VOLUME: | | RESUSCITATOR MASS: | |
|------------------------|-----------|------------------------------|-----------|
| One Hand | Two Hands | (Patient valve and bag only) | |
| Adult | 900 ml | Adult | 365 grams |
| Small Adult | 500 ml | Small Adult | 250 grams |
| Child | 360 ml | Child | 210 grams |
| Infant | 175 ml | Infant | 165 grams |

The correct ventilation frequency may vary. Please follow the current ventilation frequency recommended by the AHA.

BACKWARD AND FORWARD LEAKAGE: Negligible.
EXPIRATORY RESISTANCE: 2.1 cm H₂O @ 50 LPM
INSPIRATORY RESISTANCE: 2.9 cm H₂O @ 50 LPM

| MAXIMUM RATE (nominal): | | EXTERNAL RESUSCITATOR DIMENSIONS: | |
|---|--|-----------------------------------|--|
| Adult | 92 BPM | (Patient valve and bag only) | |
| Small Adult <td>110 BPM <td>31 cm long x 13 cm diameter</td> <td></td> </td> | 110 BPM <td>31 cm long x 13 cm diameter</td> <td></td> | 31 cm long x 13 cm diameter | |
| Child <td>168 BPM <td>32 cm long x 11.5 cm diameter</td> <td></td> </td> | 168 BPM <td>32 cm long x 11.5 cm diameter</td> <td></td> | 32 cm long x 11.5 cm diameter | |
| Infant <td>180+ BPM <td>Child 25 cm long x 9 cm diameter</td> <td></td> </td> | 180+ BPM <td>Child 25 cm long x 9 cm diameter</td> <td></td> | Child 25 cm long x 9 cm diameter | |
| | | Infant 24 cm long x 7 cm diameter | |

DEAD SPACE: 7 ml
IMPACT TEST: The Mercury Medical® CPR Bag showed no signs of damage, or degradation of performance as a result of a one meter drop onto a concrete surface.

OXIGEN TUBING CONNECTION: Universal adapter fits C.G.A. D.I.S.S. 1240 or 250 inch hose barb.
OPERATING ENVIRONMENTAL LIMITS: -18° to +50°C
STORAGE ENVIRONMENTAL LIMITS: -40° to +60°C
The CPR Bag is manufactured, approved, sponsored and endorsed by MERCURY MEDICAL® and has no connection with the Laerdal Medical Corporation.

MercuryMedical **EC/REP** **CE**
Scanlan Group B.V.
Postbus 75664
Schiphol-Tripport 1118 ZS
The Netherlands

Rx **Latex Free** **Non Sterile**

Do Not Reuse
By Prescription Only

11300 - 49th Street North
Clearwater, Florida 33762-4807
www.mercurymed.com
Patient Valve
Patent # US 5,357,951 B1

©Registered U.S. Trademark of Mercury Enterprises, Inc.
Additional copies of the "CPR Bag Directions for Use" are available on request. Contact the Marketing Dept. at 800-237-6418
Manufactured in Malaysia 04/2009 #83-900-0309 Rev. 0

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©Gedeporteerd V.S. Handelsmerk van Mercury Enterprises, Inc.