Mercury Medical
 11300 - 49th Street North
 Clearwater, Florida 33762-4807 USA
 www.mercury-med.com

Hyperinflation System (en)

with Optional Manometer/
 Pressure Relief Valve

Mercury Medical
 11300 - 49th Street North
 Clearwater, Florida 33762-4807 USA
 www.mercury-med.com

Hyperinflationssystem (da)

med optionelt manometer/
 overtryksventil

Mercury Medical
 11300 - 49th Street North
 Clearwater, Florida 33762-4807 USA
 www.mercury-med.com

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Mercury Medical
 11300 - 49th Street North
 Clearwater, Florida 33762-4807 - USA
 www.mercury-med.com

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INDICATIONS FOR USE

The Mercury Medical Hyperinflation Bag System is a pulmonary-assist device intended to provide controlled or assisted ventilation to patients.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

- This product must only be used by medical personnel well trained in techniques of positive pressure ventilation with a hyperinflation bag system.
- Never wait to begin positive pressure ventilation if a hyperinflation system is not immediately available or cannot be used effectively. (See your department manual for accepted positive pressure ventilation procedures.)
- While using supplemental oxygen, do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable chemicals.
- Should not be used in toxic or hazardous atmospheres.
- Verify proper function of the hyperinflation bag system and verify clear patient airway by monitoring for the following:
 - That the patient is being ventilated, as indicated by rise and fall of chest.
 - That the blue knurled valve rotates and the mask and patient valve are free from obstruction.
 - To clear valve obstructions, shake briskly.
 - Product should be constantly monitored while in use.
 - Gas source must be attached and flowing for hyperinflation bag system to operate properly.
 - That the breathing bag is securely and properly attached to the valve body.
- This product should be used in conjunction with a manometer pressure gauge to determine appropriate delivery pressure. Pressure may be monitored with Mercury Medical Single-Patient-Use Manometer Pressure Gauge that may be incorporated on the patient valve. Units that do not feature the built-in manometer MUST be monitored with a line attached from the pressure port to a manometer pressure gauge.
- When using the hyperinflation system, it is recommended that clinicians do not exceed 10 liters of flow per minute or 40 cm of H₂O pressure unless the user determines that higher pressures are required to effectively ventilate the patient based on clinical experience and judgment.
- Some hyperinflation units are equipped with a pressure-limiting device which opens at approximately 40 cm H₂O. These units are equipped with a lock clip. Placing the lock clip on the pressure relief valve will override the pressure limiting device. Also, an abrupt, high-volume inspiratory delivery may cause the device to exceed this level.
- Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
- Dispose of this device should be handled in accordance with hospital protocol.
- This device is assembled with PVC components which include low levels of a phthalate as indicated. Consideration should be used when treating children or pregnant or nursing women.

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DIRECTIONS FOR USE

- If a ventilation bag is not included, attach a ventilation bag to the patient valve.
- Oxygen supply tubing is attached to the oxygen inlet port. Attach opposite end of oxygen tubing to a flow regulated oxygen source.
- For units without the Mercury Medical manometer, attach pressure tubing to the pressure port and to a manometer pressure gauge. If the unit has a built-in single-patient-use manometer, pressure tubing is not necessary.
- Set the blue knurled valve to the full open position (counterclockwise).
- Adjust the oxygen flow from the gas source to the desired level while ensuring that oxygen flow is adequate for ventilatory requirements.
- Attach the patient connection to the appropriate device for delivering gas to the patient.
- Squeeze the breathing bag to deliver a breath. Release the breathing bag to allow it to refill. Observe the chest rise and fall to confirm patient is being ventilated. As bag inflates, initiate ventilation in accordance with recognized procedures.
- Adjust the blue knurled valve and oxygen flow rate as needed to provide adequate ventilation.
- Some hyperinflation units may feature a patient valve with a special pressure-limiting device mounted on the upper valve housing. If inspiration meets with pulmonary resistance, venting will occur, limiting pressure to 40cm +/-10, thereby reducing the risk of stomach distention. A hissing sound can be heard as the device opens.

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NOTE: When higher inspiration pressures are necessary, the pressure limiting device can be closed with the tip of the index finger or by engaging the lock clip while squeezing the bag.

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Mercury Medical Hyperinflation Systems Flow vs. Pressure Chart Hyperinflation Systems 0.5 and 1.0 Liter Bags			
Knob Fully Open at: Flow in Liters/Minute	(Nominal) cm of H ₂ O	Knob Fully Closed at: Flow in Liters/Minute	(Nominal) cm of H ₂ O
2 LPM	0 cm of H ₂ O	2 LPM	6.5 cm of H ₂ O
4 LPM	0 cm of H ₂ O	4 LPM	17 cm of H ₂ O
6 LPM	0.5 cm of H ₂ O	6 LPM	30 cm of H ₂ O
8 LPM	0.5 cm of H ₂ O	8 LPM	42 cm of H ₂ O
10 LPM	1 cm of H ₂ O	10 LPM	57 cm of H ₂ O
12 LPM	1.5 cm of H ₂ O	12 LPM	77 cm of H ₂ O

Mercury Medical hyperinflationssystemer, flow-tryk-tablet Hyperinflationssystemer med 0,5- og 1,0-liter-poser			
Knop helt åben ved: Flow i liter pr. minut	(Nominelt) cm H ₂ O	Knop helt lukket ved: Flow i liter pr. minut	(Nominelt) cm H ₂ O
2 lpm	0 cm H ₂ O	2 lpm	6,5 cm H ₂ O
4 lpm	0 cm H ₂ O	4 lpm	17 cm H ₂ O
6 lpm	0,5 cm H ₂ O	6 lpm	30 cm H ₂ O
8 lpm	0,5 cm H ₂ O	8 lpm	42 cm H ₂ O
10 lpm	1 cm H ₂ O	10 lpm	57 cm H ₂ O
12 lpm	1,5 cm H ₂ O	12 lpm	77 cm H ₂ O

Mercury Medical Hyperinflatiesystemen Strooming vs. druktabel Hyperinflatiesystemen zakken van 0,5 en 1,0 liter			
Knop helemaal open bij: Stroom in liter/minuut	(Nominale) cm H ₂ O	Knop helemaal dicht bij: Stroom in liter/minuut	(Nominale) cm H ₂ O
2 lpm	0 cm H ₂ O	2 lpm	6,5 cm H ₂ O
4 lpm	0 cm H ₂ O	4 lpm	17 cm H ₂ O
6 lpm	0,5 cm H ₂ O	6 lpm	30 cm H ₂ O
8 lpm	0,5 cm H ₂ O	8 lpm	42 cm H ₂ O
10 lpm	1 cm H ₂ O	10 lpm	57 cm H ₂ O
12 lpm	1,5 cm H ₂ O	12 lpm	77 cm H ₂ O

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NOTE: Nominal value only, actual output may vary depending on device, accuracy of flowmeter and the existence of tubing connection leaks.

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