

Carts/Anesthesia

Anesthesia Cart - Harloff offers a broad array of carts and accessories perfect for the demands of the anesthesia department.

MR-SAFE AND OPTIMAL



Features and Benefits

- All-aluminum Construction is lighter weight, resists organic contamination, and will not rust.
- Powder Coat Finish
- Flexible Design for Custom Drawer and Accessory Configurations
- Specially sourced ball-bearing drawer slides with low ferrous content.
- Wrap-around bumper.
- Pull out Stainless Steel side shelf.
- Plastic top with molded-in push handles
- BEST® lock for high security.
- Electronic “hands-free” locking models available.
- MR-Safe carts have been tested and certified by an independent authority.
- Delivered with wheels installed.

CLASSIC LINE



Features and Benefits

- Steel Construction with Powder-coat finish for combination of value pricing and long life.
- Ball-bearing full extension drawer slides give long life and ease of opening and closing.
- Full Vinyl wraparound bumper on separate protective steel frame.
- Carts feature standard key lock or mechanical push-button combination lock.
- Flexible Design for Custom Drawer and Accessory Configurations.
- Delivered with wheels installed.

MINI-LINE




Features and Benefits

- Space-Saving design for facilities with limited storage space.
- Steel Construction with Powder-coat finish for combination of value pricing and long life.
- Flexible Design for Custom Drawer and Accessory Configurations.
- Ball-bearing full extension drawer slides give long life and ease of opening and closing.
- Carts feature standard key lock or mechanical push-button combination lock.

Note: Product numbers indicated on pages following.

12 Powder Coat Finishes Standard



Carts/Anesthesia

Optimal™ Electronic Access Anesthesia Workstation



#GJ-PLN7

Electronic Access 7 drawer, Standard Package



#GJ-PLN7-ANS

Electronic Access 7 drawer, Specialty Package

Features and Benefits

Restricting access and key control are major issues anesthesia departments with JCAHO compliance. Now Harloff combines modern keyless security technology with functional and attractive storage cart.

- Meet or exceed all current JCAHO recommendations for the storage and transport of Anesthesia department supplies.
- Advanced Programmable Electronic Proximity Locking system uses cards, card tags, or numerical entry codes to give multiple choices for access. To gain access, the user either presents their card or enter a PIN number.
- Up to 1,000 transactions can be stored in the unit. An optional printer can print transaction log including date, time, and user.
- Keyed master lock can over-ride the electronic lock for emergency access.
- On board rechargeable power supply with recharge warning light and built-in battery saver.
- All-aluminum construction for lighter weight and durability.
- Side storage cabinet with flexible storage system for supplies or manuals.
- Full extension ball-bearing drawer slides.
- Full wraparound bumper on extended stabilizing base.
- Pull out stainless steel shelf.
- Plastic replaceable top with integrated push handles.
- 5" Sealed ball bearing casters, one directional and one braking offered.
- 31" of vertical drawer space with 2 standard drawer configurations offered. Custom drawer configurations are available.
- Full range of drawer sizes, drawer organizers, and cart accessories to meet unique department requirements.

Dimensions and Specifications

#GJ-PLN6

Drawer Sizes

1-4" (10.2 cm)
2-3.25" (8.3 cm)
2-6.5" (16.5 cm)
1-9.75" (24.8 cm)

Cart Dimensions

Standard	45.75" H	30.5" W	23.75" D
Metric	116.2 H	77.5 W	60.3 D

#GJ-PLN7

1-4" (10.2 cm)
4-3.25" (8.3 cm)
1-6.5" (16.5 cm)
1-9.75" (24.8 cm)

Net Weight

126 lbs
55.8 kg

#GJ-PLN6-ANS #GJ-PLN7-ANS

Drawer Sizes

1-4" (10.2 cm)	1-4" (10.2 cm)
2-3.25" (8.3 cm)	4-3.25" (8.3 cm)
2-6.5" (16.5 cm)	1-6.5" (16.5 cm)
1-9.75" (24.8 cm)	1-9.75" (24.8 cm)

Cart Dimensions

Standard	45.75" H	32" W	25" D
Metric	116.2 H	81.3 W	63.3 D
With Back Rail	68.5"/174 cm H		

Net Weight

139 lbs
63.1 kg

Accessories Included:

- Raised back rail system - 3 tier
- Side rails (3)
- Hospital grade 6 socket outlet with 15' cord and cord wrap
- (2) **GJ-68530-P1** drawer trays

Carts/Anesthesia/MRI

Magnetic Resonance Safe Anesthesia Carts



#GJ-MR7K-MAN

MR Safe Seven Drawer Anesthesia Cart, Specialty Package



#GJ-MR6K

MR Safe Six Drawer Anesthesia Cart, Standard Package

Features and Benefits

- Non-magnetic aluminum construction.
- Contemporary design with rounded corners.
- Specially sourced, low ferrous, ball-bearing full extension drawer slides.
- BEST® brand key lock with two keys.
- Full wrap-around vinyl bumper.
- 5" MR Safe casters, all full swivel, 2 with brake.
- Stainless steel pull-out side shelf on right side of cart.
- Replaceable plastic top with integrated push handles.
- Specific MR-Safe accessories and drawer organizers available.
- Powder coat finish - glossy white cabinet with 12 drawer colors available.
- Pre-threaded holes for easy MR Safe accessory installation.
- Cart addresses needs of hospitals and imaging centers that perform an increased variety of procedures in MRI suites.
- Provides transportation for MR Safe anesthesia machines and equipment.
- Cart is tested safe for use in both 1.5 and 3.0 Tesla magnetic resonance environments (reports are available upon request).
- Cart's durable aluminum construction resists organic contamination, will not rust, and is easily cleaned.
- Glossy white cabinet and MRI decal identifies cart as belonging in MRI suites.
- Cart is delivered with wheels installed; accessories require user installation.

Dimensions and Specifications

#GJ-MR7K-MAN #GJ-MR6K-MAN

Drawer Sizes

1-4" (10.2 cm)	1-4" (10.2 cm)
4-3.25" (8.3 cm)	2-3.25" (8.3 cm)
1-6.5" (16.5 cm)	2-6.5" (16.5 cm)
1-9.75" (24.8 cm)	1-9.75" (24.8 cm)

MR Safe Accessories Included:

- Raised back rail system - 3 tier
- Side rails (4)
- Full size drawer trays with adjustable dividers (2)

Cart Dimensions

	Standard	Metric	With Back Rail	Net Weight
Height	45.75" H	116.2 H	68.5"/174 cm H	139 lbs
Width	32" W	81.3 W		63.1 kg
Depth	25" D	63.5 D		

#GJ-MR6K

Drawer Sizes

1-4" (10.2 cm)
2-3.25" (8.3 cm)
2-6.5" (16.5 cm)
1-9.75" (24.8 cm)

#GJ-MR7K

1-4" (10.2 cm)
4-3.25" (8.3 cm)
1-6.5" (16.5 cm)
1-9.75" (24.8 cm)

Cart Dimensions

	Standard	Metric	Net Weight
Height	45.75" H	116.2 H	126 lbs
Width	30.5" W	77.5 W	55.8 kg
Depth	23.75" D	60.3 D	

Carts/Anesthesia

Optimal™ All-Aluminum Anesthesia Carts



#GJ-OPT6K

OPTIMAL Six Drawer Anesthesia Cart, Standard Package



#GJ-OPT6K-ANS

OPTIMAL Six Drawer Anesthesia Cart, Specialty Package

Features and Benefits

- Aluminum construction.
- Contemporary design with rounded corners.
- Stainless steel pull-out side shelf.
- BEST® brand key lock with two keys.
- Full wrap-around vinyl bumper.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Replaceable plastic top with integrated push handles.
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; accessories require user installation.
- Specific accessories and drawer organizers available.
- Powder coat finish - 12 colors available.
- Pre-threaded holes for easy accessory installation.
- All-aluminum construction reduces weight addressing ergonomic concerns.
- Integrated push handles on plastic top allow for additional storage on side rails.
- Top drawer is four inches high for oversized items.

Dimensions and Specifications

#GJ-OPT6K

Drawer Sizes

1-4" (10.2 cm)
2-3.25" (8.3 cm)
2-6.5" (16.5 cm)
1-9.75" (24.8 cm)

Cart Dimensions

Standard 45.75" H 30.5" W 23.75" D
Metric 116.2 H 77.5 W 60.3 D

#GJ-OPT7K

1-4" (10.2 cm)
4-3.25" (8.3 cm)
1-6.5" (16.5 cm)
1-9.75" (24.8 cm)

Net Weight

126 lbs
55.8 kg

#GJ-OPT6K-ANS

Drawer Sizes

1-4" (10.2 cm)
2-3.25" (8.3 cm)
2-6.5" (16.5 cm)
1-9.75" (24.8 cm)

Cart Dimensions

Standard 45.75" H 32" W 25" D
Metric 116.2 H 81.3 W 63.5 D
With Back Rail 68.5"/174 cm H

#GJ-OPT7K-ANS

1-4" (10.2 cm)
4-3.25" (8.3 cm)
1-6.5" (16.5 cm)
1-9.75" (24.8 cm)

Net Weight

139 lbs
63.1 kg

MR Safe Accessories Included:

- Raised back rail system - 3 tier
- Side rails (4)
- Full size drawer trays with adjustable dividers (2)
- Hospital grade 6 socket outlet with 15' cord wrap

Carts/Anesthesia

Classic Line - Malignant Hyperthermia Cart



#GJ-6640

Malignant Hyperthermia Cart, Standard Package



#GJ-6641

Malignant Hyperthermia Cart, Specialty Package

Features and Benefits

- Some patients under anesthesia experience a sudden and dangerous rise in blood temperature, a condition called Malignant Hyperthermia (MH). Healthcare facilities that administer volatile inhalational agents known to trigger MH must have emergency equipment immediately available for the treatment of MH (Accreditation Association for Ambulatory Care). Surgery centers must have immediately available, and in an organized cart, all items required for treatment of MH (American Society of Anesthesiologists and the American Association of Nurse Anesthetists).
- The Harloff Malignant Hyperthermia Cart will demonstrate JCAHO compliance and/or best practice for all OR suites, Level 1 Trauma Units, Surgery Centers, and Ambulatory Facilities.
- Cart includes: integrated 1.6 cu ft refrigerator/freezer with hospital grade plug, 10 glycolpacks, LED illuminated digital thermometer with cart power safety feature, 3" and 6" storage drawers, stainless steel breakaway locking system on both the refrigerator and the two storage drawers, full wrap-around vinyl bumper on separate steel frame, ball-bearing full extension drawer slides, removable plastic top, 5-inch premium wheels with one brake and one tracking caster, powder coat finish in 12 colors.
- The cart increases patient safety and departmental efficiency by providing a convenient, portable and economic solution for storage and transport of cooled glycol packs, cold saline IV solutions, and dantrolene (Dantrium) used to stabilize MH. Cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most likely needed.
- External digital thermometer is visible at a distance in all lighting conditions, and allows continual monitoring of the supplies' temperature.
- Cart ships with wheels assembled, so the cart does not need to be inverted to install the wheels, thereby protecting the refrigeration components.

Dimensions and Specifications

#GJ-6640

Drawer Sizes

1-3"

1-6"

2-Compartment with pre-installed refrigerator

Cart Dimensions Net Weight

Standard 41.5" H 32" W 22" D 146 lbs

Metric 105.4 H 81.3 W 55.9 D 66.4 kg

#GJ-6641

Drawer Sizes

1-3"

1-3"

2-Compartment with refrigerator

Cart Dimensions

Standard 44.5" H 32" W 22" D 154 lbs

Metric 113.0 H 81.3 W 45.7 D 70 kg

Accessories Included:

- IV Pole (extends to 85" / 216 cm)
- Full Drawer Tray Organizer
- Side-Mounted Drop Shelf
- Stainless Steel top Rail/Handle

Carts/Anesthesia

Classic Line - Accessorized Workstation Carts



#GJ-6456

Deluxe Six Drawer Workstation Cart with key lock



#GJ-7456

Deluxe Six Drawer Workstation Cart with mechanical combination pushbutton lock

Features and Benefits

- Key lock with two keys (#GJ-6456)
- Mechanical combination pushbutton lock (#GJ-7456).
- Hospital grade 6 socket electrical outlet with 15' cord and cord wrap.
- Full wrap-around vinyl bumper on separate steel frame.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Removable plastic top.
- Raised back rail system with 2 horizontal rails.
- Aluminum push handle/side rail (4).
- Two 3-inch exchange trays and two 6 inch exchange trays.
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; some accessories require user installation.
- Full range of accessories and drawer organizers available.
- Unibody steel construction.
- Powder coat finish - 12 colors available.
- Combination lock models eliminate the need for keys; combinations can be reset by user.
- Cart is accessorized for simplified ordering.
- Cabinet and drawer design allows for customized drawer configurations.

Dimensions and Specifications

#GJ-6456 Accessories shown in photo are included

Drawer Sizes

4-3" (7.6 cm)

1-6" (15.2 cm)

1-2" (30.5 cm)

Cart Dimensions

Standard 41.5" H 32" W 22" D 171 lbs

Metric 105.4 H 81.3 W 55.9 D 77.6 kg

with back rail 57.75"/146.7 cm H

#GJ-7456 Accessories shown in photo are included

Drawer Sizes

4-3" (7.6 cm)

1-6" (15.2 cm)

1-12" (30.5 cm)

Cart Dimensions

Standard 42.5" H 32" W 22" D 171 lbs

Metric 108 H 81.3 W 55.9 D 77.6 kg

with back rail 58.25"/148 cm H

Carts/Anesthesia

Classic Line - Tall Workstation Carts, Combination Lock



#GJ-7451

Six Drawer Workstation Cart with mechanical combination pushbutton lock, Specialty Package



#GJ-7450

Six Drawer Workstation Cart with mechanical combination pushbutton lock, Standard Package

Features and Benefits

- Mechanical combination pushbutton lock (**GJ-7451**).
- Full wrap-around vinyl bumper on separate steel frame.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Removable plastic top.
- Raised back rail system with 2 horizontal rails.
- Aluminum push handle/side rail (4).
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; some accessories require user installation.
- Full range of accessories and drawer organizers available.
- Unibody steel construction.
- Powder coat finish - 12 colors available.
- Combination lock models eliminate the need for keys; combinations can be reset by user.
- Pre-installed wheels reduce set-up time when delivered.
- Cabinet and drawer design allows for customized drawer configurations.
- Drawers are double-walled with no sharp edges.

Dimensions and Specifications

#GJ-7451

Drawer Sizes

4-3" (7.6 cm)
1-6" (15.2 cm)
1-12" (30.5 cm)

Cart Dimensions

Standard	42.5" H	32" W	22" D
Metric	108 H	81.3 W	55.9 D
with back rail	58.25"/148 cm H		

Net Weight

140 lbs
63.5 kg

#GJ-7450

Drawer Sizes

4-3" (7.6 cm)
1-6" (15.2 cm)
1-12" (30.5 cm)

Cart Dimensions

Standard	42.5" H	32" W	22" D
Metric	108 H	81.3 W	55.9 D

Net Weight

145 lbs
65.8 kg

Carts/Anesthesia

Classic Line - Tall Workstation Carts



#GJ-6450

Six Drawer Workstation Cart with key lock,
Standard Package



#GJ-6451

Six Drawer Workstation Cart with key lock,
Specialty Package

Features and Benefits

- Key lock with two keys.
- Full wrap-around vinyl bumper on separate steel frame.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Removable plastic top.
- Aluminum push handle/side rail (2).
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; some accessories require user installation.
- Full range of accessories and drawer organizers available.
- Unibody steel construction.
- Powdercoat finish - 12 colors available.
- Combination lock models eliminate the need for keys; combinations can be reset by user.
- Pre-installed wheels reduce set-up time when delivered.
- Cabinet and drawer design allows for customized drawer configurations.
- Drawers are double-walled with no sharp edges.

Dimensions and Specifications

#GJ-6450

Drawer Sizes

- 4-3" (7.6 cm)
- 1-6" (15.2 cm)
- 1-12" (30.5 cm)

Cart Dimensions

	41.5" H	32" W	22" D	Net Weight
Standard	41.5" H	32" W	22" D	145 lbs
Metric	105.4 H	81.3 W	55.9 D	65.8 kg

#GJ-6451

Drawer Sizes

- 4-3" (7.6 cm)
- 1-6" (15.2 cm)
- 1-12" (30.5 cm)

Cart Dimensions

	41.5" H	32" W	22" D	Net Weight
Standard	41.5" H	32" W	22" D	140 lbs
Metric	105.4 H	81.3 W	55.9 D	63.6 kg
with back rail	57.75"/146.7 cm H			

Carts/Anesthesia

Classic Line - Short Workstation Carts, Combination Lock



#GJ-7351

Four Drawer Workstation Cart with combination pushbutton lock, Specialty Package



#GJ-7350

Four Drawer Workstation Cart with mechanical combination pushbutton lock, Standard Package

Features and Benefits

- Mechanical combination pushbutton lock.
- Full wrap-around vinyl bumper on separate steel frame.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Removable plastic top.
- Raised back rail system with 2 horizontal rails.
- Aluminum push handle/side rail (4).
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; some accessories require user installation.
- Full range of accessories and drawer organizers available.
- Unibody steel construction.
- Powder coat finish - 12 colors available.
- Combination lock models eliminate the need for keys; combinations can be reset by user.
- Pre-installed wheels reduce set-up time when delivered.
- Cabinet and drawer design allows for customized drawer configurations.
- Drawers are double-walled with no sharp edges.

Dimensions and Specifications

#GJ-7351

Drawer Sizes

- 2-3" (7.6 cm)
- 1-6" (15.2 cm)
- 1-12" (30.5 cm)

Cart Dimensions

	Standard	Metric	with back rail	Net Weight
Height	36.5" H	92.7 H	58.25"/148 cm H	115 lbs
Width	32" W	81.3 W		55.9 D
Depth	22" D	55.9 D		52.2 kg

#GJ-7350

Drawer Sizes

- 2-3" (7.6 cm)
- 1-6" (15.2 cm)
- 1-12" (30.5 cm)

Cart Dimensions

	Standard	Metric	Net Weight
Height	36.5" H	92.7 H	120 lbs
Width	32" W	81.3 W	55.9 D
Depth	22" D	55.9 D	54.5 kg

Carts/Anesthesia

Classic Line - Short Workstation Carts



#GJ-6350

Four Drawer Workstation Cart with key lock,
Standard Package



#GJ-6351

Four Drawer Workstation Cart with key lock,
Specialty Package

Features and Benefits

- Key lock with two keys.
- Full wrap-around vinyl bumper on separate steel frame.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Removable plastic top.
- Aluminum push handle/side rail (2).
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; some accessories require user installation.
- Full range of accessories and drawer organizers available.
- Unibody steel construction.
- Powder coat finish - 12 colors available.
- Combination lock models eliminate the need for keys; combinations can be reset by user.
- Pre-installed wheels reduce set-up time when delivered.
- Cabinet and drawer design allows for customized drawer configurations.
- Drawers are double-walled with no sharp edges.

Dimensions and Specifications

#GJ-6350

Drawer Sizes

2-3" (7.6 cm)
1-6" (15.2 cm)
1-12" (30.5 cm)

Cart Dimensions

Standard	35.5" H	32" W	22" D	Net Weight	120 lbs
Metric	90.2 H	81.3 W	55.9 D		54.5 kg

#GJ-6351

Drawer Sizes

2-3" (7.6 cm)
1-6" (15.2 cm)
1-12" (30.5 cm)

Cart Dimensions

Standard	35.5" H	32" W	22" D	Net Weight	115 lbs
Metric	90.2 H	81.3 W	55.9 D		52.2 kg
with back rail	51.5"/130.8 cm H				

Carts/Anesthesia

Classic Line - Monitor Carts



#GJ-6650

with optional GJ-680404 electrical outlet and cord wrap
Six drawer Monitor Cart with key lock, Standard package



#GJ-6651

Six Drawer Monitor Cart with
key lock, Specialty Package

Features and Benefits

- Short overhead shelf system with 2 shelves (#GJ-6650).
- Tall overhead shelf system with 2 shelves offering multiple mounting locations (#GJ-6651).
- Key lock with two keys.
- Full wrap-around vinyl bumper on extended width steel frame.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Removable plastic top.
- Aluminum push handle/side rail (1).
- Hospital grade 6 socket electrical outlet with 15' cord and cord wrap (#GJ-6651).
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; some accessories require user installation.
- Full range of accessories and drawer organizers available.
- Unibody steel construction.
- Powder coat finish - 12 colors available.
- Extended bumper provides greater stability when shelves are loaded.
- Cabinet and drawer design allows for customized drawer configurations.
- Drawers are double-walled with no sharp edges.

Dimensions and Specifications

#GJ-6650

Drawer Sizes	Cart Dimensions	Net Weight
3-3" (7.6 cm)	Standard 41.5" H 40" W 25" D	180 lbs
2-6" (15.2 cm)	Metric 105.4 H 63.5 W 55.9 D	81.7 kg
1-9" (22.9 cm)	with Monitor Shelves 66.5"/168.9 cm H	

#GJ-6651

Drawer Sizes	Cart Dimensions	Net Weight
3-3" (7.6 cm)	Standard 41.5" H 40" W 25" D	185 lbs
2-6" (15.2 cm)	Metric 105.4 H 101.6 W 63.5 D	84 kg
1-9" (30.5 cm)	with Monitor Shelves 75"/190.5 cm H	

Carts/Anesthesia

Classic Line - Secure Storage Cart



#GJ-6610

Eight Drawer Storage Cart with key lock,
Standard Package



#GJ-6611

Eight Drawer Storage Cart with key lock,
Specialty Package

Features and Benefits

- Key lock with two keys.
- Full wrap-around vinyl bumper on separate steel frame.
- 5" Sealed ball-bearing premium casters, all full swivel, 2 with brake.
- Aluminum push handle/side rail (#GJ-6610/#GJ-6600).
- Stainless Steel top rail/push handle (#GJ-6611/#GJ-6601).
- Removable stainless steel top (#GJ-6601).
- Removable plastic top.
- Ball-bearing full extension drawer slides.
- Cart is delivered with wheels installed; accessories require user installation.
- Full range of accessories and drawer organizers available.
- Unibody steel construction.
- Powdercoat finish - 12 colors available.
- Pre-installed wheels reduce set-up time when delivered.
- Cabinet and drawer design allows for customized drawer configurations.
- Drawers are double-walled with no sharp edges.

Dimensions and Specifications

#GJ-6610

Drawer Sizes
6-3" (7.6 cm)
2-6" (15.2 cm)

#GJ-6611

Drawer Sizes
6-3" (7.6 cm)
2-6" (15.2 cm)

Accessories Included:

- Utility Hooks (4)
- Stainless steel top rail
- 3 - 6 Compartment trays
- 3 - 7 Compartment trays

#GJ-6600

- Same specifications as **#GJ-6610** with break-away lock.

#GJ-6601

- Same specifications as **#GJ-6611**.

Cart Dimensions	Net Weight
Standard 41.5" H 32" W 22" D	155 lbs
Metric 105.4 H 81.3 W 55.9 D	70.4 kg

Cart Dimensions	Net Weight
Standard 44.5" H 32" W 22" D	125 lbs
Metric 113.0 H 81.3 W 55.9 D	56.8 kg

Carts/Anesthesia

Mini Line - Workstation Carts



#GJ-3145K-ANS

Five Drawer Short Mini-Line Workstation Cart with key lock, Specialty package



#GJ-3156K-ANS

Six Drawer Tall Mini Line Workstation Cart with key lock, Specialty Package

Features and Benefits

- Narrow design with reduced footprint.
- Two cabinet heights - with vertical drawer space of 24" (small) or 30" (tall).
- Key lock with two keys.
- Ball-bearing full extension drawer slides.
- All steel construction.
- 3" Stem casters, all full swivel, 2 with brake
- Removable plastic top.
- Powder coat finish - 12 colors available.
- Specific accessories and drawer organizers available.
- Space-saver design makes cart ideal for anesthesiologists, clinics, doctor offices, phlebotomists, and special procedures.
- Two sizes accommodate a variety of storage requirements.
- Narrow cart design allows supply storage in areas where space is limited.

Dimensions and Specifications

#GJ-3145K-ANS

Drawer Sizes	Cart Dimensions	Net Weight	Accessories Included:
3-3" (7.6 cm)	Standard 32" H 18" W 18" D	88 lbs	<ul style="list-style-type: none"> • Raised back rail system - 2 tier • Side rails (2) • 6 Tilt bin organizer
1-6" (15.2 cm)	Metric 81.3 H 45.7 W 45.7 D	40.0 kg	
1-9" (22.9 cm)	with Monitor Shelves 45.63"/123.5 cm H		

#GJ-3156K-ANS

Drawer Sizes	Cart Dimensions	Net Weight	Accessories Included:
3-3" (7.6 cm)	Standard 41.5" H 32" W 22" D	108 lbs	<ul style="list-style-type: none"> • Raised back rail system - 2 tier • Side rails (2) • 6 Tilt bin organizer
2-6" (15.2 cm)	Metric 105.4 H 81.3 W 55.9 D	49 kg	
1-9" (30.5 cm)	with Monitor Shelves 54.75"/139.1 cm H		

Mini Line - Workstation Carts



#GJ-3145K
Five Drawer Short Mini-Line Workstation Cart
with key lock, Standard package



#GJ-3156K
Six Drawer Tall Mini Line
Workstation Cart with
key lock, Standard Package

Features and Benefits

- Narrow design with reduced footprint.
- Two cabinet heights - with vertical drawer space of 24" (short) or 30" (tall).
- Key lock with two keys.
- Ball-bearing full extension drawer slides.
- All steel construction.
- 3" Stem casters, all full swivel, 2 with brake
- Removable plastic top.
- Powder coat finish - 12 colors available.
- Specific accessories and drawer organizers available.
- Space-saver design makes cart ideal for anesthesiologists, clinics, doctor offices, phlebotomists, and special procedures.
- Two sizes accommodate a variety of storage requirements.
- Narrow cart design allows supply storage in areas where space is limited.

Dimensions and Specifications

#GJ-3145K

Drawer Sizes	Cart Dimensions	Net Weight
3-3" (7.6 cm)	Standard 32" H 18" W 18" D	75 lbs
1-6" (15.2 cm)	Metric 81.3 H 45.7 W 45.7 D	34.1 kg
1-9" (22.9 cm)		

#GJ-3156K

Drawer Sizes	Cart Dimensions	Net Weight
3-3" (7.6 cm)	Standard 38.25" H 18" W 18" D	95 lbs
2-6" (15.2 cm)	Metric 97.2 H 45.7 W 45.7 D	43.2 kg
1-9" (30.5 cm)		

Carts/Anesthesia

E-Series Line



#GJ-E30-6K
Standard 6 drawer 60" cabinet



#GJ-E24-4K
Standard 4 drawer 24" cabinet

Features and Benefits

- Value priced line with Harloff quality.
- Key lock with two keys.
- Drawer fronts have integrated pull handle for easy opening anywhere on the face.
- Ball Bearing Drawer Slides for smooth opening over the life of the cart.
- Corner bumpers to protect walls.
- Removable plastic top.
- Aluminum push handle.
- 5" casters, all full swivel, one with brake, one directional for better steering control.
- Cart is delivered with wheels installed to save time putting into service.
- Powder Coat finish with 12 colors available.
- Full range of accessories and drawer organizers available.

Dimensions and Specifications

#GJ-E30-6K

Drawer Sizes	Cart Dimensions	Net Weight
4-3" (7.6 cm)	Standard 41.5" H 32" W 22" D	145 lbs
1-6" (15.2 cm)	Metric 105.4 H 81.3 W 55.9 D	65.8 kg
1-12" (30.5 cm)		

#GJ-E24-4K

Drawer Sizes	Cart Dimensions	Net Weight
2-3" (7.6 cm)	Standard 35.5" H 32" W 22" D	120 lbs
1-6" (15.2 cm)	Metric 90.2 H 81.3 W 55.9 D	54.5 kg
1-12" (30.5 cm)		


Carts/Accessories

Magnetic Resonance Safe




#GJ-MR-RBR3

Raised Back Rail System - 3 Tier for mounting accessories.



#GJ-MR-RBR2

Raised Back Rail System - 2 Tier for mounting accessories.



#GJ-MR-SR

Side rail for mounting accessories.



#GJ-MR-RLHK

Rail Hook, fits all accessory mounting rails.



#GJ-MR-4TILT 4-Compartment Tilt Bin

#GJ-MR-5TILT 5-Compartment Tilt Bin


#GJ-MR-6TILT 6-Compartment Tilt Bin

All include Mounting Clips



#GJ-MR-02

Stainless Steel oxygen tank holder for D/E size cylinders



#GJ-CARDBRD

Cardiac board with 2 brackets



#GJ-MR-IV2

2 prong IV pole, adjustable height



#GJ-68530-P1 3" Full drawer tray with adjustable dividers

Fit one per drawer

#GJ-40130 3W Drawer divider tray only

#GJ-68530-P2 3" Full drawer tray with extra adjustable dividers



#GJ-680503

9" and 12" drawer divider set with 3 long and 2 short dividers
Dims: 17"L x 22.75"W x 8.75"H
Fit one per drawer



#GJ-40109 6-compartment Drawer Tray

Dims: 14.5"L x 7.6"W x 1.75H

Fit three per drawer

#GJ-40108 6 compartment drawer tray

Dims: 14.5"L x 7.6"W x 1.75"H

Fit three per drawer

Optimal Line



#GJ-OPT-RBR3

Raised Back Rail System - 3 Tier for mounting accessories.



#GJ-OPT-RBR2

Raised Back Rail System - 2 Tier for mounting accessories.



#GJ-OPT-SR

Side Rail for mounting accessories.



#GJ-OPT-RLHK

Rail Hook, fits all accessory mounting rails.



#GJ-OPT-OUTLET

Hospital Grade 6-socket Electrical Outlet with 15' cord and cord wrap.



#GJ-AL2362

Glove box holder, aluminum



#GJ-OPT-POSHELF

Stainless Steel pull-out shelf, replacement only



#GJ-OPT-02

Stainless Steel oxygen tank holder for D/E size cylinders

Carts/Accessories

Classic Line - E-Series



#GJ-680311
Short Aluminum Overhead Monitor Support Set - (shelves sold separately)

#GJ-680311FI

Short Aluminum Overhead Set Support Set For Field Installation

#GJ-680312

14" Overhead Monitor Shelf

#GJ-680313

20" Overhead Monitor Shelf

#GJ-680314

Stainless Steel 14" Overhead Monitor Shelf



#GJ-680310
Tall Aluminum Overhead Monitor Support Set

(shelves sold separately)

#GJ-680310FI

Tall Aluminum Overhead Support Set - For Field Installation

#GJ-680317

Rubber Mat



#GJ-68032
Adjustable overhead Shelf with Accessory Bar



#GJ-68163
Raised Back Rail System-2 Rails

#GJ-68157

Raised Back Rail System-3 Rails



#GJ-680414

Utility Box


#GJ-680424

Utility Box with Rail Clip



#GJ-680419

Rail Clip



#GJ-680431

Small Utility Basket with Nylon Mounting Clamps

#GJ-40256

Small Utility Basket (13"L x 4"D x 4.5"H)



#GJ-680447

Large Utility Basket with Mounting Brackets

#GJ-680407

Large Utility Basket with Rail Clips

#GJ-40257

Large Utility Basket only (10"L x 7.5"D x 5.25"H)



#GJ-680415

Tape/Label Dispenser, includes two rail clips for mounting on push handle or raised back rail system



#GJ-680413

18" Catheter Box with rail clip, for mounting on push handle

#GJ-680434

36" Catheter Box with rail clip



#GJ-683803

4-Compartment Tilt Bin


#GJ-683804

5-Compartment Tilt Bin

#GJ-683805

6-Compartment Tilt Bin

All Include Mounting Clips



#GJ-680404

Hospital Grade 6-Socket electrical Outlet, with 15' cord & cord wrap

#GJ-680405

Cord Wrap Only



#GJ-680428

Three gallon Waste Container and mounting bracket

#GJ-40717

Three Gallon Waste Container only

#GJ-40638

Replacement mounting bracket for three gallon waste container



#GJ-680408

Two Gallon Plastic Waste Container includes mounting hardware



#GJ-680409

Two Gallon Metal Waste Container with Cover includes mounting hardware

#GJ-680410

Two Gallon Metal Waste Container without Cover Includes mounting hardware



#GJ-68161

Aluminum Push Handle/Side Rail, 18"

#GJ-68162

Aluminum Back Rail, 24" for Classic and E-Series Only

Carts/Accessories

Classic Line - E-Series Drawer Organizers



#GJ-681501
3" Adjustable Exchange Tray with
3 Long and 3 Short Dividers
Dims: 16.5"L x 11"W x 2.25"H
Fit two per drawer

#GJ-40179
3" Tray Long Dividers, Pk 3

#GJ-40180
3" Tray Short Dividers, Pk-3



#GJ-681502
6" Adjustable Exchange Tray with
3 Long and 3 Short Dividers
Dims: 16.5"L x 11"W x 5.25"H
Fit two per drawer

#GJ-40181
6" Tray Short Dividers, Pk 3

#GJ-40182
6" Tray Long Dividers, Pk-3



#GJ-68530-P1
3" Full Drawer Tray
with Adjustable Dividers
Fit one per drawer

#GJ-40582
3" Full Drawer Tray only

#GJ-68530-P2
3" Full Drawer
tray with extra Adjustable Dividers



#GJ-68605
5 Space Divider
Insert pack of 10 Fits all 3"
Full Drawer trays

#GJ-68604
4 Space Divider Insert
Pack of 10

#GJ-68603
3 Space Divider Insert
Pack of 10

#GJ-68602
2 Space Divider
Insert pack of 10



#GJ-40345

Utility Bin

#GJ-40361

Utility Bin Dividers, pack of 10



#GJ-680503
9" and 12" Drawer Divider
set with 3 long
and 2 short dividers
Dims: 17"L x 22.75"W x 8.75"H
Fit one per drawer



#GJ-40109
7-Compartment Drawer Tray
Dims: 14.5"L x 7.6"W x 1.74H
Fit three per drawer

#GJ-40108
6 Compartment Drawer Tray
Dims: 14.5"L x 7.6"W x 1.75"H
Fit three per drawer

Mini Line



#GJ-MINI-RBR
Raised Back Rail System - For
MINI Line cart



#GJ-683803
4-Compartment Tilt Bin
#GJ-683804
5-Compartment Tilt Bin
#GJ-683805
6-Compartment Tilt Bin
All Include Mounting Clips



#GJ-680404
Hospital Grade 6-Socket
electrical Outlet, with 15'
cord & cord wrap
#GJ-680405
Cord Wrap Only



#GJ-68161
Aluminum Push
Handle/Side Rail, 18'



#GJ-68620-P1
Mini Line Drawer Tray
with Adjustable Dividers
Fit one per drawer
#GJ-68620-P2
Mini Line Drawer Tray
with extra Adjustable Dividers



#GJ-OPT-RLHK
Rail Hook, fits all
accessory mounting
rails.