



# Endobronchial Tubes

## Rusch Double Lumen PVC Endobronchial Tubes

The very best endobronchial tubes available.


Rusch tubes are supplied with secure high quality connections, wide port openings, and clear color coding throughout.

Product	Description	Item Number	Size/I.D.	Qty
Robert Shaw For Left Bronchus Sterile	• Ruschelit® PVC	#12-11610026	26F	1
	• Single-use	#12-11610028	28F	1
	• X-ray opaque stripe	#12-11610035	35F	1
	• Graduated	#12-11610037	37F	1
	• Stylet included	#12-11610039	39F	1
	• Color coded for easy identification of ports, balloons and 15mm connectors	#12-11610041	41F	1



Carlens Tube Sterile	• Ruschelit® PVC	#12-11610135	35F	1	For Left Bronchus
	• Single-use	#12-11610137	37F	1	
	• X-ray opaque stripe	#12-11610139	39F	1	
	• Graduated	#12-11610141	41F	1	
	• Stylet included				

- Color coded for easy identification of ports, balloons and 15mm connectors
- Low pressure/high volume cuffs
- Components as shown in diagram
- With carina hook



Robert Shaw For Right Bronchus Sterile	• Ruschelit® PVC	#12-11620026	26F	1	For Right Bronchus
	• Single-use	#12-11620028	28F	1	
	• X-ray opaque stripe	#12-11620035	35F	1	
	• Graduated	#12-11620037	37F	1	
	• Stylet included	#12-11620039	39F	1	
	• Color coded for easy identification of ports, balloons and 15mm connectors	#12-11620041	41F	1	

- Low pressure/high volume cuffs
- Components a shown in diagram



Dr. White For Right Bronchus Sterile	• Ruschelit® PVC	#12-116201350	35F	1	For Right Bronchus
	• Single-use	#12-116201370	37F	1	
	• X-ray opaque stripe	#12-116201390	39F	1	
	• Graduated	#12-116201410	41F	1	
	• Stylet included				

- Color coded for easy identification of ports, balloons and 15mm connectors
- Low pressure/high volume cuffs
- Components a shown in diagram
- With carina hook

