

Masks

MERCURY MEDICAL® CLEAR ANATOMICAL FACE MASKS (Child, Pediatric, Infant, Neonatal)

The anatomically-designed contoured cushion conforms to natural facial contours allowing a secure low pressure seal and reduces damage to the eyes. Special attention has been given to the face cone of the mask raising the contour over the bridge of the nose. This accommodates pronounced facial features and eliminates contact with the patient. Clear construction allows continuous visual patient contact and reduces patient anxiety.

Features:

- Completely disposable. Eliminates cross-contamination risks and associated reprocessing costs.
- Clear material construction for better patient visualization.
- Anatomically designed. Reduces patient discomfort and accommodates a wide range of facial features.
- Soft, pliable, air-filled cushion. Reduces wrinkles or creases and allows for a secure, low pressure seal.
- Removable color-coded headstrap hookrings. Allows use of traditional headstraps. Removable for use in hand-held procedures. (Infant and Neonatal masks supplied without head straps)




SPECIFICATIONS

Child -	Accommodates most child needs.**
Pediatric -	Accommodates most toddlers and small children.**
Infant -	Accommodates most infants.***
Neonatal -	Accommodates most newborns.***

** Masks have 22mm Connecting Port (meets ISO Standards)

***Masks have 15mm Connecting Port (meets ISO Standards)



#10-50500	Adult Large Clear Anatomical Face Mask	(Large)	30/box
#10-50501	Adult Small Clear Anatomical Face Mask	(Medium)	30/box
#10-50502	Child Clear Anatomical Face Mask	(Small)	20/box
#10-50503	Pediatric Clear Anatomical Face Mask		20/box
#10-50504	Infant Clear Anatomical Face Mask		20/box
#10-50505	Neonatal Clear Anatomical Face Mask		20/box

Masks

Face Masks


Pediatric

Features anatomically correct "average" contour for infant faces. Provides overall tight seal with internal cavity, minimizing dead space. Designed specifically for children under the age of 10. 22mm female fitting

RBS STYLE

Plastic - non-conductive

#12-111300	Size 0 (premature)	2/Box
#12-111301	Size 1 (up to 1 year)	2/Box
#12-111302	Size 2 (1-3 years)	2/Box
#12-111303	Size 3 (4-10 years)	2/Box




LATEX
FREE

RBS STYLE CONDUCTIVE

Rubber - non-conductive

#12-310300	Size 0 (premature)
#12-310301	Size 1 (up to 1 year)
#12-310302	Size 2 (1-3 years)
#12-310303	Size 3 (4-10 years)*

*size 3 mask supplied with hook ring




Contains
LATEX

ROUND TYPE SILICONE MASKS

15mm male fitting for size 0 and 1,
22mm fitting for size 2


#10-50601	Size 0 (48mm)	10/pkg
#10-50602	Size 1 (60mm)	10/pkg
#10-50603	Size 2 (76mm)	10/pkg



LATEX
FREE

ANATOMICAL MASKS

D #07-860220	Child Silicone Mask No. 3-4 w/Multi-Function Mask Cover
E #07-870220	Adult Silicone Mask No. 4-5+ w/Multi-Function Mask Cover
F #07-865200	Multi-Function Mask Cover for Child Mask 3-4
G #07-875200	Multi-Function Mask Cover for Adult Mask 4-5+



Masks

AMBU TYPE

Non-conductive transparent masks ideal for resuscitation, obstetrics, emergencies and anesthesia where conductivity is not essential, but where transparency is advantageous.

#19-E1720	Infant #0
#19-E1720A	Infant #0A 15mm male
#19-E1722	Infant #2
#19-E1724	Infant #4
#19-E1725	Infant #5



#19-E1720A




#19-E1725

ANATOMICAL FACE MASK

A face mask designed in an anatomical shape. The medical air cushion inflating valve facilitates a tight seal. Available in sizes from (very small) infant to (extra large) adult. A five prong hook ring is provided.

#12-3101-100	Infant	#0
#12-3101-101	Small Child	#1
#12-3101-102	Child	#2
#12-3101-103	Small Adult	#3
#12-3101-104	Medium Adult	#4
#12-3101-105	Adult	#5
#12-3101-106	Large Adult	#6



#06-769122
Plugs For Mask

CRYSTALONE™ FACE MASK

Reusable, lightweight mask for resuscitation and anesthesia. The strong, crystal-clear polysuphone cone allows for early regurgitation or aspiration and for the anesthetist to monitor lip color and intermittent misting. The super soft inflatable pad is specifically designed for a comfortable, safe seal.

#12-3105102	Child/Small Adult (medium) w/hook ring
#12-3105103	Adult (large) w/hook ring




#12-3105102

Masks


ECO Mask™ Clear Non-PVC Anesthetic Face Masks

Eco Mask™ is a range of environmentally friendly anesthetic face masks that are anatomically shaped so that they fit the contours of the patient's face giving a precise seal. Eco Mask utilizes a soft silicone-like seal that gives improved patient comfort and the bonded transparent shell provides excellent visual inspection during use. There are seven sizes available in the range.




Extra Large Adult Anesthetic Face Mask Size 6

#DM-7096 | Box Quantity: 50




Large Adult Anesthetic Face Mask Size 5

#DM-7095 | Box Quantity: 50




Medium Adult Anesthetic Face Mask Size 4

#DM-7094 | Box Quantity: 50




Small Adult Anesthetic Face Mask Size 3

#DM-7093 | Box Quantity: 50



Pediatric Anesthetic Face Mask Size 2

#DM-7092 | Box Quantity: 50



Infant Anesthetic Mask Size 1

#DM-7091 | Box Quantity: 50



Neonatal Anesthetic Face Mask Size 0

#DM-7090 | Box Quantity: 50



Mask Sprays/Head Harnesses

Maskumm™ Anesthesia Mask Spray - Economically Scents Masks For Easier Pediatric Induction.

- Maskumm™ eliminates the need for costly manufactured prescented masks.
- Facilitates mask acceptance by pediatric patients.
- Two popular scents, cherry and bubblegum.
- Parent and child become actively involved in the induction by selecting the cherry or bubblegum.
- Fast and easy to apply.
- Maskuum™ does not require refrigeration.
- Environmentally safe, does use fluorocarbons.
- Many adults prefer Maskumm™ with approximately 300 sprays per bottle.



Item Number **#56-250BG** Bubblegum scent
Item Number **#56-250C** Cherry scent

HEAD HARNESES

MASK HARNESS

- Molded, not stamped
- Reinforced holes on adult & child sizes for durability
- Anti-static, latex free neoprene
- Reusable
- Center vents for comfort and fit at back of head



Item No	Description
#53-4-1054-21	Infant Head Harness
#53-4-1054-22	Child Head Harness
#53-4-1054-23	Adult Head Harness