

SpO₂ Adapter - and Extension cables

Quick Reference Guide and Connector Chart
































































































Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Connector	Connector to Monitor Description	Connector
AMI	all (Nellcor based)	FQ-92/CA032320	8 / 2,4	DB9 female		7pin male Hypertronics	
Aquitron	Serena MX Sleep System (BCI based)	FQ-92/CA013320	8 / 2,4	DB9 female		7pin male Redel	
Armeda	PM 9000 Express, BX66-2808Y	FQ-92/CA841320	8 / 2,4	DB9 female		6pin male Redel	
Artema/S+W	Gamme, Athena, Diascope 2 for use with Nellcor	FQ-92/CA071320	8 / 2,4	DB9 female		10pin female Nicolay	
Athena	Diascope 2, Gamme 8000 (Nellcor SpO ₂)	FQ-92/CA074320	8 / 2,4	DB9 female		10pin female Nicolay	
Baxter/Colin	9050/9060	FQ-92/CA075320	10 / 3,0	On request, please call sales office!			
	ASAT to Nellcor	FQ-92/CA081320	8 / 2,4	DB9 female		9pin male Hypertronics	
	Explorer, Monarch	FQ-92/CE082320	8 / 2,4	DB9 female		DB9 male	
BCI (Biochem)	3040, 3100, 3101, 3301, 3303 3401, 9100, 9200 Advisor, Autocorr, FingerPrint, Mini Torr, Mini Torr Plus, 6100 only (Serena)	FQ-92/CA091320	8 / 2,4	DB9 female		7pin male Redel	
	PM 4002	FQ-92/CA093120	4 / 1,2	DB9 female		5pin male Redel	
	3180, 3300, 3301, 3302, 3303, 3304, 3401/-2/-3, CapnoCheck/II/sleep/plus	FQ-92/CE092320	10 / 3,0	DB9 female		DB9 male	
Bird	4400 (P/N CSI 511-SD)	FQ-92/CA111320	8 / 2,4	DB9 female		5pin male Lemo	
Bionet	BM31	FQ-92/CE852320	8 / 2,4	DB9 female		7pin bionet male	
Bruker (Odam)	all Odam	FQ-92/CA121320	8 / 2,4	DB9 female		12pin male Nicolay	
Burdick (Spacelabs)	100, PRO2	FQ-92/CE131320	8 / 2,4	DB9 female		DB9 male	
Colin	BP-88/306/508/808/885, Press-mate SAT, Advantage 9861/9871, PressMate-PAL 3118, Prodigy II 2210/ 2210-P/2240/2240-P (Nellcor based)	FQ-92/CE141320	10 / 3,0	DB9 female		DB9 male	
Corometrics (GE Marquette)	188, 128, 128F, 129, 129F, 511, 556 (Nellcor based)	FQ-92/CA151320	8 / 2,4	DB9 female		12pin male Hypertronics	
Corpuls	Life Pack, 08/16, (BCI based)	FQ-92/CA161320	8 / 2,4	DB9 female		12pin male Nicolay	
	Corpuls monitors	FQ-92/CA162320	10 / 3,0	DB9 female		12pin male Nicolay	
Criticare (CSI)	503, 504, 507, 508, 581, 1100, Poet, Poni, Scholar	FQ-92/CA171320	10 / 3,0	DB9 female		5pin male Lemo	
	503 Spot, 503DX, 504DX, 506, 506DXNP, mini SpO ₂	FQ-92/CE172320	10 / 3,0	DB9 female		DB9 male	
Critikon/ Dinamap	Dinamap 9700, 9710, 9720 plus, Compact. MPS, Select	FQ-92/CE182320	10 / 3,0	DB9 female		DB9 male	
	Dinamap, 8700, 8710, 8720, Plus, Oxyshuttle 1/2	FQ-92/CA181320	10 / 3,0	6pin Jack Compushield		7pin Lemo	
	Dinamap, 8700, 8710, 8720, Plus, Oxyshuttle 1/2	FQ-92/CE183320	8 / 2,4	6pin Jack Compushield		6pin Compushield plug	
	Oxyshuttle 1/2	FQ-92/CA184320	10 / 3,0	6pin Jack Compushield	12pin female Thorkom		
	Dinamap Pro300/400, Pro 1000	FQ-92/CA185320	10 / 3,0	DB9 female		14 pin plug 3M	
Datascope	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascope SpO ₂), Passport (Datascope SpO ₂),	FQ-92/CA192320	8 / 2,4	DB9 female		8 pin male DIN	
	Expert, Passport (Nellcor SpO ₂)	FQ-92/CE191320	10 / 3,0	DB9 female		DB9 male	
	Accutorr Plus (Nellcor SpO ₂)	FQ-92/CA193320	8 / 2,4	DB9 female		14 pin plug 3M	

Chart continued on following page...

SpO₂ Adapter - and Extension cables

Quick Reference Guide and Connector Chart

...Chart continued from previous page.






















































Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Connector	Connector to Monitor Description	Connector
Datex Ohmeda	AS-3; CS-3; CH/S, CH/2S, CH-RS, CG/2S, Cardiocap I, II, 5; Capnomac Ultima; Light Monitor; M-Estpr, OscarII; Satellite Trans, Plus, Plus, ULTSVI, ULTSVO, ULTSIO ULTS00	FQ-92/CA204320	10 / 3,0	DB9 female		10 pin female Nicolay	
	3700, 3710, 3770, 3775, 3800, 3900, Biox3740, 3760, 4700, OxiCap, 5250RGM, Rascal II, TuffSat	FQ-92/CA201320	8 / 2,4	DB9 female		9 pin male Hypertronics	
Digicare	Digimax 5500 (BCI SpO2)	FQ-92/CA231320	8 / 2,4	DB9 female		8pin male DIN	
Dinamap/ Critikon	please see Critikon, GE						
Drager	Vitalert 1000, Cicero	FQ-92/CA253320	8 / 2,4	DB9 female		8pin male Lemo	
	Oxipac 2000, 2500 (BCI based)	FQ-92/CE251320	10 / 3,0	DB9 female		DB9 male	
	Narkomed	FQ-92/CA252320	8 / 2,4	DB9 female		12 pin male Lemo	
Entri Medical	Propaq 102, 104, 106 (Nellcor based)	FQ-92/CA771320	8 / 2,4	DB9 female		12 pin male Lemo	
GE Marquette	DASH, Eagle 3000, 4000, Solar, TRAM X50-Modules (Nellcor SpO2)	FQ-92/CA286320	10 / 3,0	DB9 female		11 pin male blue flat	
	Nellcor Oxysmart XL	FQ-92/CA282320	10 / 3,0	DB9 female		11 pin male blue flat	
	7000 Series TRAM A.R., TRAM, Modular-X00SL Series (Ohmeda based)	FQ-92/CA281320	10 / 3,0	DB9 female		9 pin male Hypertronics	
Hellige	Eagle 1000 (Nellcor based)	FQ-92/CA285320	10 / 3,0	DB9 female		10pin male	
Goldway	UT 4000 A (BCI based)	FQ-92/CA791320	8 / 2,4	DB9 female		5 pin male Redel	
	UT 4000 A (BCI based)	FQ-92/CA792120	4 / 1,2	DB9 female		5 pin male Redel	
Healthdyne	All (Nonin SpO2)	FQ-92/CE311320	10 / 3,0	DB9 female		DB9 male	
Hewlett Packard	please see Philips Medical						
Innomed		FQ-92/CA331320	8 / 2,4	DB9 female		7 pin male Redel	
Invivo Research	4500 Plus 2 (Nellcor SpO2)	FQ-92/CA351320	8 / 2,4	DB9 female		7 pin male Lemo	
IVY Biomedical	405P, 405T	FQ-92/CE361320	8 / 2,4	DB9 female		DB9 male	
	405A, 405D	FQ-92/CA363320	8 / 2,4	DB9 female		12pin male Lemo	
Kontron	Micromon, Minimon, Supermon, Kolormon, 7840, 7845, 7137	FQ-92/CA371320	8 / 2,4	DB9 female		12 / 6 pin male Nicolay	
Lynn Medical	3300, 3301, 3302, 3303, 3304, 3401, (BCI SpO2)	FQ-92/CE381320	8 / 2,4	DB9 female		DB9 male	
Magna	Mon-Mtek-30, 33, 35, 40, 80, 85, 1200, 1201, 1300, 1401, 1406, 1411, 1413, 1501, 1503, MOD02, MOD03 (BCI SpO2)	FQ-92/CE391320	8 / 2,4	DB9 female		DB9 male	
Med Air	LifeSense, PulsSense (BCI SpO2)	FQ-92/CE431320	8 / 2,4	DB9 female		DB9 male	
Mediana	YM-2002 Series (Nellcor SpO2)	FQ-92/CA801320	8 / 2,4	DB9 female		14 pin plug 3M	
M.D.E.	Escort II OPT30, Prism.	FQ-92/CE422320	8 / 2,4	DB9 female		DB9 male	
	Escort 100-OPT 11A, Escort 300A	FQ-92/CA421320	8 / 2,4	DB9 female		12pin male Lemo	

...Chart continued on following page.

SpO₂ Adapter - and Extension cables

Quick Reference Guide and Connector Chart

...Chart continued from previous page.

























































Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Connector	Connector to Monitor Description	Connector
Medlab	POX 10 u. POX 10L monitor elan, /C, VITRO	FQ-92/CA442220	5,3 / 1,6	DB9 female		7pin male Redel	
	POX 10 u. POX 10L monitor elan, /C, VITRO	FQ-92/CE443220	5,3 / 1,6	7pin female Redel		7pin male Redel	
	NANOX2 (BCI SpO2)	FQ-92/CE441320	5,3 / 1,6	DB9 female		DB9 male	
Medsource	BPM 200 (Nellcor SpO2)	FQ-92/CE811320	8 / 2,4	DB9 female		DB9 male	
MEK	all disposables at MP-100, -400, -1000	FQ-92/CE821320	8 / 2,4	DB9 female		DB9 male	
Mennen	Horizon XL, MR1300, Mercury (Nellcor SpO2)	FQ-92/CA451320	8 / 2,4	DB9 female		10pin male Lemo	
	Envoy	FQ-92/CA452320	8 / 2,4	DB9 female			
Mindray	PM 7000, 8000, 9000	FQ-92/CA461320	8 / 2,4	DB9 female		6pin male Redel	
	PM 8000 Express	FQ-92/CA462320	8 / 2,4	DB9 female		6pin male Redel	
MSA/Catalyst Research	Mini OX IV	FQ-92/CA481320	8 / 2,4	DB9 female		7pin male Lemo	
Nellcor/Mallinckrodt	N-100 (w/High Sat), N-200 (pre-amp)	FQ-92/CA491320	8 / 2,4	DB9 female		12pin male Lemo	
	NPB-290/5, N-390, N-3000, N6000 (w/Oxismart)	FQ-92/CA496320	10 / 3,0	DB9 female		14pin plug 3M	
	for use in conjunction with OxiMax Sensors	FQ-92/CA498220	6,7 / 2,0	DB9 female		14pin plug 3M	
	All (not listed above)	FQ-92/CE497220	4 / 1,2	DB9 female		DB9 male	
Nihon Kohden	All (not listed above)	FQ-92/CE497320	10 / 3,0	DB9 female		DB9 male	
	AL-800PA Module; BSM-1101, 2101/2, 3101, 7100, 7200; OLV-1100A, Lifescope	FQ-92/CA501320	8 / 2,4	DB9 female		10pin male Nicolay	
	AL900P Modul; BSM 1102, 9510, OGS 2000 (Nellcor SpO2)	FQ-92/CE502320	10 / 3,0	DB9 female		DB9 male	
	BSM-4101, 4103, 4111, 4113 Life scope P; BSM-5105, 5135 Life scope A; BSM-2351, 2353 Life scope L; BSM-2301, 2303 Life scope I	FQ-92/CA504320	8 / 2,4	9 pin female Nihon Kohden		14 pin male Nihon Kohden	
	BSM-4104A/ 4114A (Nellcor SpO2)	FQ-92/CA503320	8 / 2,4	DB9 female		14pin plug 3M	
Nonin	Avant 2120, 9600, 9700, 8500 Series, 9480 Series	FQ-92/CE512320	8 / 2,4	DB9 female		DB9 male	
	8500, 8600, 8700, 8800, 9600, 9700	FQ-92/CA513320	8 / 2,4	DB9 female		AMP 8pin	
	8604 Oxygen / Biochem 1040A	FQ-92/CA511320	8 / 2,4	DB9 female		6pin male Lemo	
Novamatrix/Respironics	505, 510, 515B/C, 520, 7000, 7100, Oxypleth superbright	FQ-92/CE521320	8 / 2,4	7pin female Hypertronics		7pin male Hypertronics	
	505, 510, 515B/C, 520, 7000, 7100, Oxypleth superbright	FQ-92/CA522320	8 / 2,4	DB9 female		7pin male Hypertronics	
Odam	please see Bruker						
Ohmeda	3700, 3710, Biox 3740/3760, 3770, 3775, 3800, 4700 OxiCap, Rascal II	FQ-92/CA541320	8 / 2,4	DB9 female		9 pin male Hypertronics	
Pace Tech	Vitalmax, 800 series (Nellcor SpO2)	FQ-92/CA571320	8 / 2,4	DB9 female		5pin female Switchcraft	
Penlon	PM 9000 (Nellcor SpO2)	FQ-92/CA831320	8 / 2,4	DB9 female		6pin male Redel	

...Chart continued on following page.

SpO₂ Adapter - and Extension cables

Quick Reference Guide and Connector Chart

...Chart continued from previous page.

Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Connector Connector	Connector to Monitor Description	Connector Connector
Philips Medical	78352A/C, 78354A/C, 78834C, Merlin: 1020A, 1024A, 1025A/B, 1094A, 1205A, 1350B, 1722A/B, 1732A/B, 2475B, Miracle system	FQ-92/CA323320	10 / 3,0	DB9 female		12pin male Merlin	
	Philips IntelliVue M3001A, Module M 1020B, Heart Stream/HeartStart XL	FQ-92/CA322320	10 / 3,0	DB9 female		8pin male Viridia	
	Philips IntelliVue M3001A, Module M 1020B, Heart Stream/HeartStart XL	FQ-92/CA322120	3,6 / 1,1	DB9 female		8pin male Viridia	
	same as above (Nellcor SpO2)	FQ-92/CA321220	6,7 / 2,0	8pin female Viridia		12pin male Merlin	
	Philips A1: M3921A, M2923A Philips A3: M3927A, M3929A	FQ-92/CE324320	10 / 3,0	DB9 female		DB9 male	
	HP Viridia, M3, M4; Philips IntelliVue M3001 series	FQ-92/CE325220	6,7 / 2,0	8pin female Viridia		8pin male Viridia	
Physio Control	Lifepak 20	FQ-92/CA591320	8 / 2,4	DB9 female		14pin plug 3M	
	Lifepak 12 (Nellcor SpO2)	FQ-92/CE592320	10 / 3,0	DB9 female		DB9 male	
Protocol	Propaq 206-EL, 204-EL, 202-EL, 106-EL (Nellcor SpO2)	FQ-92/CE601320	8 / 2,4	DB9 female		DB9 male	
RGB Medical	Omicrom FT-/PREMIUM (BCI SpO2)	FQ-92/CA611320	8 / 2,4	DB9 female		7pin male Redel	
Schiller	Oximeter OX-1 (Nellcor SpO2)	FQ-92/CA623320	8 / 2,4	DB9 female		14pin plug 3M	
	Argus Pro, TM-7 (Nellcor SpO2)	FQ-92/CA622320	8 / 2,4	DB9 female		7pin male DIN	
	CM 8 (BCI SpO2)	FQ-92/CA621320	8 / 2,4	DB9 female		6pin male DIN	
	Physiogard TM 910, (BCI based)	FQ-92/CA624320	8 / 2,4	DB9 female		12pin male Nicolay	
Sensormedics	Oxyshuttle I, II	FQ-92/CA631320	10 / 3,0	6pin jack Compushield		12pin female Thorkom	
	Oxyshuttle I, II	FQ-92/CE632320	8 / 2,4	6pin jack Compushield		6pin plug Compushield	
		FQ-92/CA633320	10 / 3,0	6pin jack Compushield		7pin male Lemo	
Siemens	MicrO2: SC 8000, 9000; Sirecust 404-1, 700, 630; Sirem SpO2 Module (Nellcor SpO2)	FQ-92/CA643320	8 / 2,4	DB9 female		10pin male Nicolay	
	Sirecust 9000, SC 6000, 7000	FQ-92/CA641120	4 / 1,2	DB9 female		7pin male Siemens	
	Sirecust 9000, SC 6000, 7000	FQ-92/CA642220	7 / 2,2	DB9 female		7pin male Siemens	
Spacelabs	90419-01, 90651A-07	FQ-92/CA667320	8 / 2,4	DB9 female		7pin male Hypertronics	
	90419-01, 90501, 90651A-07, 90351-D/Gm; 90465; 90466 (Novamatrix SpO2)	FQ-92/CE668320	8 / 2,4	7pin female Hypertronics		7pin male Hypertronics	
	90351-0/6, 90465, 90466, 90467, 90489, 90651A, IM77; PC Express (Nellcor SpO2)	FQ-92/CA661320	8 / 2,4	DB9 female		12pin male Hypertronics	
		FQ-92/CA662320	8 / 2,4	DB9 female		12pin male Hypertronics	
	90651-08, 90496 Ultraview (Nellcor)	FQ-92/CA665320	8 / 2,4	DB9 female		10pin female Spacelabs	
	PC Scout, 90303B (Nellcor based)	FQ-92/CA663320	8 / 2,4	DB9 female		12pin male Hypertronics	
	90651-08, 90496 Ultraview (Novamatrix)	FQ-92/CA666020	0,67 / 0,2	7pin female Hypertronics		10pin female Spacelabs	
	90419-01, 90501, 90651A-07, 90351-D/Gm; 90465; 90466 (Novamatrix SpO2)	FQ-92/CA664020	0,67 / 0,2	7pin female Hypertronics		12pin male Hypertronics	


For any manufacturer not mentioned above please contact us. No responsibility is taken for the correctness of this information!
All manufacturers' trademarks and registered names recognized. This chart is subject to change without prior notice.

Note for understanding the abbreviation:
CA = Adapter Cable
CE = Extension Cable

SpO₂ Reusable Fingerprobes

Quick Reference Guide and Connector Chart

Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Picture
Aequitron	Serena MX 9101-OX	FQ-92011320	10 / 3	7 pin male Redel	1
		FQ-92012320	10 / 3	9 pin male Hypertronics	2
Alaris	Vital Check 4400 (Criticare based)	FQ-92731120	3 / 0.9	DB9 male	36
	Vital Check 4400 (Criticare based)	FQ-92732320	10 / 3	DB9 male	36
Allegiance Healthcare	Allegiance	FQ-92021320	10 / 3	14 pin 3M plug	4
		FQ-92022120	3 / 0.9	DB9 male	3
Alliance Instruments	Integriti Model	FQ-92741120	3 / 0.9	7 pin male Hypertronics	5
AMI	all (Nelcor based)	FQ-92032120	3 / 0.9	DB9 male (must use adapter)	3
	all	FQ-92031320	10 / 3	7 pin male Hypertronics	5
AMS	MOM Toco, MOM US	FQ-92751320	10 / 3	5 pin male Lemo	6
Anamed	All	FQ-92041320	10 / 3	7 pin male Redel	1
Aquitron	Serena MX Sleep System (BCI based)	FQ-92761120	3 / 0.9	DB9 male (must use adapter)	3
	Serena MX Sleep System 9101-OX (Ohmeda based)	FQ-92762320	10 / 3	7 pin male Redel	1
		FQ-92763320	10 / 3	9 pin male Hypertronics	2
Armeda	PM 9000 Express	FQ-92841320	10 / 3	6 pin male Redel	37
	PM 9000 Express	FQ-92842120	3 / 0.9	DB9 male (must use adapter)	3
Armstrong	AD900/1000/2000/3000 (BCI based)	FQ-92061120	3 / 0.9	DB9 male	32
Artema / S+W Athena	MM200	FQ-92071320	10 / 3	12 pin male Nicolay	7
	Diascope G2, Gamme 8000 (Nelcor based)	FQ-92072120	3 / 0.9	DB9 male (must use adapter)	3
	S&W 9050/9060	FQ-92073320	10 / 3	10 pin female Nicolay	8
	Diascope G2, Gamme 8000 (Nelcor based)	FQ-92075320	10 / 3	10 pin female Nicolay	35
Atom Medical	V505/900/2001/Atom 2000	FQ-92051320	10 / 3	14 pin 3M plug	4
Baxter/Colin	Asat	FQ-92081320	10 / 3	9 pin male Hypertronics	2
	Explorer/Monarch (Nelcor based)	FQ-92082120	3 / 0.9	DB9 male	3
BCI Smith Medical (Biochem)	3180, 3300, 3301, 3302, 3303, 3304, 3401/-2/-3, CapnoCheck/II/sleep/plus	FQ-92091120	3 / 0.9	DB9 male	32
	3100, 3101, 6100 Serena, 9100, 9200 Advisor, Autocorr/-Plus, FingerPrint, Mini Torr, Mini Torr Plus		3 / 0.9	DB9 male (must use adapter)	32
	3180, 3300, 3301, 3302, 3303, 3304, 3401/-2/-3, CapnoCheck/II/sleep/plus	FQ-92092320	10 / 3	DB9 male	32
	3100, 3101, 6100 Serena, 9100, 9200 Advisor, Autocorr/-Plus, FingerPrint, Mini Torr, Mini Torr Plus	FQ-92093320	10 / 3	7 pin male Redel	1
Biosys	BPM200 (Nelcor based)	FQ-92101120	3 / 0.9	DB9 male	3
Bird	4400 (P/N CSI 511-10L)	FQ-92111320	10 / 3	5 pin male Lemo	6
	4400 (P/N CSI 511-SD)	FQ-92112120	3 / 0.9	DB9 male	36
Bruker	4044, 3444 (BCI based)	FQ-92121120	3 / 0.9	DB9 male (must use adapter)	32
Burdick	100, PRO ₂ (BCI based)	FQ-92131120	3 / 0.9	DB9 male	32
Colin	BP-88/306/508/808/885, Press-mate SAT, Advantage 9861/9871, PressMate- PAL 3118, Prodigy II 2210/2210-P/2240/ 2240-P (Nelcor based)	FQ-92141120	3 / 0.9	DB9 male	3
	BP 308/408, BX5	FQ-92142320	10 / 3		
	BX 10ma	FQ-92143320	10 / 3	14 pin 3M plug	4
Corometrics (GE Marquette)	118, 120, 510	FQ-92151320	10 / 3	12 pin male Hypertronics	28
	511 (Ohmeda based)	FQ-92153320	10 / 3	9 pin male Hypertronics	2
	188, 128, 128F, 129, 129 F, 511, 556 (Nelcor based)	FQ-92152120	3 / 0.9	DB9 male (must use adapter)	3
Corpulse	Life Pack, 08/16 (BCI based)	FQ-92161120	3 / 0.9	DB9 male (must use adapter)	32
Criticare (CSI)	501, 503, 504, 507, 508, 581, 1100, POET (602-3A), PONI, Scholar, 8100/8500 Poet	FQ-92171320	10 / 3	5 pin male Lemo	6
	503 Spot, 503DX, 504DX, 506, 506DXNP, mini SpO ₂	FQ-92172120	3 / 0.9	DB9 male	36
	503 Spot, 503DX, 504DX, 506, 506DXNP, mini SpO ₂	FQ-92173320	10 / 3	DB9 male	36
	503 Spot, 503DX, 504DX, 506, 506DXNP, mini SpO ₂	FQ-92174320	10 / 3	DB9 male rectangular	36




...Chart continued on following page.

SpO₂ Reusable Fingerprobes

Quick Reference Guide and Connector Chart

...Chart continued from previous page.

Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Picture
Critikon/Dinamap	Dinamap Oxytrak, 8280	FQ-92181120	3 / 0.9	DB9 male	32
	Dinamap, 8700, 8710, 8720, Plus, Oxyshuttle 1/2	FQ-92182120	3 / 0.9	6 pin Compushield plug (must use adapter)	9
	Dinamap, 8700, 8710, 8720, Plus	FQ-92183320	10 / 3	7 pin male Lemo	10
	Dinamap 9700, 9710, 9720 plus, Compact. MPS, Pro300/400, Select (Nellcor based)	FQ-92184120	3 / 0.9	DB9 male	3
	Dinamap, 8700, 8710, 8720, Plus, Oxyshuttle 1/2	FQ-92185320	10 / 3	12 pin female Thorkom	39
Datascope	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascope based), Passport 2/EL2 (Datascope based)	FQ-92191320	10 / 3	8 pin male DIN	11
	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascope based), Passport (Datascope based)	FQ-92192120	3 / 0.9	DB9 male (must use adapter)	30
	Expert, Passport (Nellcor based)	FQ-92193120	3 / 0.9	DB9 male	3
	Accutorr Plus (Nellcor based), Passport XG	FQ-92195320 FQ-92196320	10 / 3 10 / 3	14 pin 3M plug 14 pin 3M plug	4 4
Datex Ohmeda see also Ohmeda	AS-3; CS-3; CH/S, CH/2S, CH-RS, CG/2S, Cardicap I, II, 5; Capnomac Ultima; Light Monitor; M-Estpr, Oscar II; OSP-200, Satellite Trans, Plus, ULTSVI, ULTSVO, ULTSIO, ULTS00	FQ-92201320	10 / 3	10 pin female Nicolay	12
	same as above	FQ-92202120	3 / 0.9	DB9 male (must use adapter)	32
	same as above	FQ-92202220	6.6 / 2	DB9 male (must use adapter)	32
	3700, 3710, 3770, 3775, 3800, 3900, Biox3740, 3760, 4700, OxiCap, 5250RGM, Rascal II, TuffSat	FQ-92203120	3 / 0.9	9 pin male Hypertronics	2
	3700, 3710, 3770, 3775, 3800, 3900, Biox3740, 3760, 4700, OxiCap, 5250RGM, Rascal II, TuffSat	FQ-92203320	10 / 3	9 pin male Hypertronics	2
Digicare	Digioxii 920, PO 930, Oxitrace, Digimax	FQ-92231320	10 / 3	8 pin male DIN	34
	Digimax 5500 (BCI based)	FQ-92233120	3 / 0.9	DB9 male	32
Dima Italy	all (Nellcor based)	FQ-92211120	3 / 0.9	DB9 male	3
Dinamap/Critikon	Dinamap, 8700, 8710, 8720, Plus, Oxyshuttle 1/2	FQ-92221120	3 / 0.9	6 pin Compushield plug	9
	Dinamap, 8700, 8710, 8720, Plus	FQ-92222320	10 / 3	7 pin male Lemo	10
	Dinamap 9700, 9710, 9720, plus, Compact. MPS, Pro300/400, Select (Nellcor based)	FQ-92223120	3 / 0.9	DB9 male	3
	Dinamap Oxytrak, 8280	FQ-92224120	3 / 0.9	DB9 male	32
	Dinamap, 8700, 8710, 8720, Plus, Oxyshuttle 1/2	FQ-92225320	10 / 3	12 pin female Thorkom	39
Dixtal	DX 2010 series, DX 8100	FQ-92241320	10 / 3	7 pin male Hypertronics	5
Drager	Oxipac 2000, 2500 (BCI based)	FQ-92251120	3 / 0.9	DB9 male	32
	Vitalert 1000, Cicero, Narkomed, Vitara, PM8014, 8060 (Nellcor based)	FQ-92252120	3 / 0.9	DB9 male	3
	Narkomed	FQ-92253320	10 / 3	12 pin male Lemo	14
Elmed	PSO1000 (BCI based)	FQ-92261120	3 / 0.9	DB9 male	32
Entri Medical	Propaq 102, 104, 106 (Nellcor based)	FQ-92771120	3 / 0.9	DB9 male	3
	Propaq 102, 104, 106	FQ-92772320	10 / 3	12 pin male Lemo	14
Fukuda Denshi	HG302 Module, DS-5100E, DS-5300W (Nellcor based)	FQ-92271120	3 / 0.9	DB9 male	3
Gambro	BP806	FQ-92291120	3 / 0.9	DB9 male	3
GE Marquette	DASH, Eagle 3000, 4000, Solar, TRAM X-50-Modules (Nellcor based)	FQ-92281120	3 / 0.9	DB9 male (must use adapter)	3
	7000 Series TRAM A.R., TRAM, Modular-X00SL Series (Ohmeda based)	FQ-92282320	10 / 3	9 pin male Hypertronics	2
	DASH, Eagle, Solar, TRAM 250, 450, 650, 850 Modules (straight through)	FQ-92283320	10 / 3	11 pin male blue flat	15
	Handheld (Nonin based)	FQ-92284120	3 / 0.9	DB9 male	32
Hellige	Eagle 1000, 2000, 3000	FQ-92285120	3 / 0.9	DB9 male	3
	SMU, SMK, SMV	FQ-92286320	10 / 3	9 pin male Hypertronics	2
Getemed	VitalGuard (Nellcor based)	FQ-92301120	3 / 0.9	DB9 male	3
	VitalGuard	FQ-92302320	10 / 3	14 pin 3M plug	4
Goldway	UT 6000 A/F (Nellcor based)	FQ-92791120	3 / 0.9	DB9 male	3
	UT 4000 A (BCI based)	FQ-92792320	10 / 3	5 pin male Redel	31
Hellige	please see GE Marquette				




...Chart continued on following page.

SpO₂ Reusable Fingerprobes

Quick Reference Guide and Connector Chart

...Chart continued from previous page.

Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Picture
Healthdyne	All (Nonin based)	FQ-92311120	3 / 0.9	DB9 male	32
Hewlett Packard/ Agilent	please see Philips Medical				
Innomed	Oxycard, Innocare-T (Nonin based)	FQ-92331120	3 / 0.9	DB9 male	32
	(Nonin based)	FQ-92332220	5 / 1.5	7 pin Redel male	41
Invacare	IRC 700, 725, 735, 775 (BCI based)	FQ-92341120	3 / 0.9	DB9 male	32
Invivo Research	4500 Scout, 4500 Plus, Scout+/-1/+2/+3/+4, Scout ES, Millennia, Omega 1445, Omni-Trak TVS, Sat Scout	FQ-92351320	10 / 3	7 pin male Lemo	10
	Same as above and 4500 Plus 2 (Nellcor based)	FQ-92352120	3 / 0.9	DB9 male (must use adapter)	3
IVY Biomedical	405A, 405D (Nellcor based)	FQ-92361120	3 / 0.9	DB9 male	3
	Sat-Guard Ivy2000, 405TM, 405CM	FQ-92362320	10 / 3	14 pin 3M plug	4
	405P, 405T	FQ-92361120	3 / 0.9	DB9 male (must use adapter)	3
Kontron	Micromon, Minimon, Supermon, Kolomon, 7840, 7845, 7137	FQ-92371320	10 / 3	12 pin male Nicolay	19
	same as above (Nellcor based)	FQ-92372120	3 / 0.9	DB9 male	3
Lynn Medical	3300, 3301, 3302, 3303, 3304, 3401, (BCI based)	FQ-92381120	3 / 0.9	DB9 male	32
Magna	Mon-Mtek-30, 33, 35, 40, 80, 85, 1200, 1201, 1300, 1401, 1406, 1411, 1413, 1501, 1503, MOD02, MOD03 (BCI based)	FQ-92391120	3 / 0.9	DB9 male	32
Mallinckrodt	please see Nellcor				
Marquette	please see GE Marquette				
M.D.E.	Escort 100-OPT 11A, Escort 300A, Escort II OPT 30, Prism. (Nellcor based)	FQ-92421120	3 / 0.9	DB9 male	3
Med Air	LifeSense, PulsSense (BCI based)	FQ-92431120	3 / 0.9	DB9 male	32
Mediana	YM-2002 Series (Nellcor based)	FQ-92801120	3 / 0.9	DB9 male (must use adapter)	3
Medlab	NANOX2 (BCI based)	FQ-92441120	3 / 0.9	DB9 male	32
	NANOX2 (BCI based)	FQ-92443320	10 / 3	DB9 male	32
	POX 10 u. POX 10L, elan, /C Monitor, VITRO	FQ-92442320	10 / 3	7 pin male Redel	40
Medsource	BPM 200 (Nellcor based)	FQ-92811120	3 / 0.9	DB9 male	3
MEK	MP-100, -110, -111, -400, -800, -1000	FQ-92821320	10 / 3	DB9 male	18
Mennen	Horizon XL (Ohmeda based)	FQ-92451320	10 / 3	9 pin male Hypertronics	2
	Horizon XL, MR1300, Mercury	FQ-92452320	10 / 3	10 pin male Lemo	38
	Horizon XL, MR1300, Mercury, Envoy (Nellcor based)	FQ-92453120	3 / 0.9	DB9 male (must use adapter)	3
Mindray	PM-509B, 600, 5000, 6000, 8000, 9000 (Nellcor based)	FQ-92461120	3 / 0.9	DB9 male	3
	PM 7000/8000/9000, PM 600	FQ-92462320	10 / 3	6 pin male Redel	37
		FQ-92463320	10 / 3	12 pin male Merlin	16
M.R.L.	All	FQ-92471320	10 / 3	5 pin female Switchcraft	20
MSA/Catalyst Research	Mini OX IV	FQ-92481320	10 / 3	7 pin male Lemo	10
	Mini OX IV	FQ-92482120	3 / 0.9	DB9 male (must use adapter)	3
	Mini OX V	FQ-92483320	10 / 3	7 pin male Lemo	10
Nellcor/ Mallinckrodt	NPB-290/5, N-390, N-3000, N6000 (w/Oxismart)	FQ-92491320	10 / 3	14 pin 3M plug	4
	N-100 (w/High Sat), N-200	FQ-92492120	3 / 0.9	DB9 (must use pre-amp)	3
	All other (not listed above) w/o Oximax technology	FQ-92493120	3 / 0.9	DB9 male	3
	All other (not listed above) w/o Oximax technology	FQ-92494320	10 / 3	DB9 male	3
Nihon Kohden	AL-800PA Module; BSM-1101, 2101/2, 3101, 7100, 7200; OLV-1100A, Lifescope	FQ-92501320	10 / 3	10 pin male Nicolay	21
	same as above	FQ-92502120	3 / 0.9	DB9 male (must use adapter)	32
	AL900P Modul, BSM 1102, 9510, OGS 2000, BSM-4104A/4114A (Nellcor based)	FQ-92503120	3 / 0.9	DB9 male	3
	BSM2301K, 2303K, 4103, 4113, 9510J/K	FQ-92505120	3 / 0.9	9 pin male Nihon	23
	BSM2301K, 2303K, 4103, 4113, 9510J/K	FQ-92504220	5 / 1.6	9 pin male Nihon	23



























...Chart continued on following page.

SpO₂ Reusable Fingerprobes

Quick Reference Guide and Connector Chart

...Chart continued from previous page.

Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Picture	
Nihon Kohden	BSM-4101, 4103, 4111, 4113 Life scope P; BSM-5105, 5135 Life scope A; BSM-2351, 2353 Life scope L; BSM-2301, 2303 Life scope I	FQ-92506320	10 / 3	14 pin male Nihon Kohden		1
Nonin	Avant 2120, 9600, 9700, 8500 Series, 9480 Series, 8600	FQ-92511120	3 / 0.9	DB9 male		3
	Avant 2120, 9600, 9700, 8500 Series, 9480 Series, 8600	FQ-92513320	10 / 3	DB9 male		3
	8604 Oxygen/Biochem 1040A	FQ-92512320	10 / 3	6 pin male Lemo		24
Novamatrix/Respironics	500 (non-superbright)	FQ-92521320	10 / 3	7 pin male Hypertronics		5
	505, 510, 515B/C, 520, 7000, 7100, Oxypleth superbright (same as above)	FQ-92522320	10 / 3	7 pin male Hypertronics		5
	Handheld 512, 513, 710, 715	FQ-92523120	3 / 0.9	DB9 male (must use adapter)		32
Respironics	920M/Plus (Nonin based)	FQ-92524120	3 / 0.9	DB9 male		25
Odam	All (BCI based)	FQ-92531120	3 / 0.9	DB9 male (must use adapter)		32
Ohmeda see also Datex-Ohmeda	2000, 405T	FQ-92541320	10 / 3	14 pin 3M plug		4
	3700, 3710, Biox 3740/3760, 3770, 3775, 3800, 4700 OxiCap, Rascal II, TuffSat	FQ-92542320	10 / 3	9 pin male Hypertronics		2
		FQ-92542120	3 / 0.9	9 pin male Hypertronics		2
		FQ-92543120	3 / 0.9	DB9 male (must use adapter)		32
Omnicon	RGB Monitor	FQ-92551320	10 / 3	7 pin male Redel		1
Ortivus	Mida Bedside	FQ-92561320	10 / 3	10 pin female Nicolay		8
Pace Tech	All / 520	FQ-92571320	10 / 3	5 pin female Switchcraft		20
	Vitalmax, 800 series (Nellcor based)	FQ-92572120	3 / 0.9	DB9 male (must use adapter)		3
Palco	All	FQ-92581120	3 / 0.9	RJ-12 plug		26
	All	FQ-92582320	10 / 3	RJ-12 plug		26
Penlon	PM 9000 (Nellcor based)	FQ-92831120	3 / 0.9	DB9 male (must use adapter)		3
	PM 9000	FQ-92832320	10 / 3	6 pin male Redel		37
Philips Medical	78352A/C, 78354A/C, 78834C, Merlin: 1020A, 1024A, 1025A/B, 1094A, 1205A, 1350B, 1722A/B, 1732A/B, 2475B, Miracle system.	FQ-92321320	10 / 3	12 pin male Merlin		16
	same as above (Nellcor based)	FQ-92322120	3 / 0.9	DB9 male (must use adapter)		3
	Philips Monitor A1/A3 (Nellcor based)			DB9 male		
	Philips Monitor A1/A3 (Nellcor based)	FQ-92322220	6,7 / 2	DB9 male		3
	HP Viridia, M3, M4, Philips M3000A Series, Philips IntelliVue M3001A, Module M1020B, Heart Stream/ HeartStart XL (Nellcor based)	FQ-92327120	3 / 0.9	8 pin male HP Viridia		17
		FQ-92327220	5 / 1,6			
	FQ-92327320	10 / 3				
Physio Control	Lifepak 12 (Nellcor based)	FQ-92591120	3 / 0.9	DB9 male		3
	Lifepak 20	FQ-92592320	10 / 3	14 pin 3M plug		4
	Lifestat 1600	FQ-92593320	10 / 3	9 pin male Hypertronics		42
Protocol	Propaq 206-EL, 204-EL, 202-EL, 106-EL (Nellcor based)	FQ-92601120	3 / 0.9	DB9 male		3
RGB Medical	Omicrom FT/-PREMIUM	FQ-92611320	10 / 3	7 pin male Redel		1
	Omicrom FT/-PREMIUM (BCI based)	FQ-92612120	3 / 0.9	DB9 male (must use adapter)		32
Schiller	Oximeter OX-1, Argus Pro	FQ-92621320	10 / 3	14 pin 3M plug		4
	Oximeter OX-1 / OX-2, Argus Pro, TM-7 (Nellcor based)	FQ-92622120	3 / 0.9	DB9 male (must use adapter)		3
	CM 8 (BCI based)	FQ-92623120	3 / 0.9	DB9 male (must use adapter)		32
Sensormedics	Oxyshuttle I, II	FQ-92631120	3 / 0.9	6 pin Compushield plug (must use adapter)		9
	Oxyshuttle I, II	FQ-92632320	10 / 3	12 pin female Thorkom		39
Siemens	MicroO2: SC 6000, 7000, 8000, 9000; Sirecust 404-1, 700, 630; Sirem SpO2 Module	FQ-92641320	10 / 3	10 pin male Siemens		27
	MicroO2: SC 6000, 7000, 8000, 9000; Sirecust 404-1, 700, 630; Sirem SpO2 Module (Nellcor based)	FQ-92642120	3 / 0.9	DB9 male (must use adapter)		3
Simed	S100, S50, S-100E	FQ-92651320	10 / 3	9 pin male Hypertronics		2
Spacelabs	90419, 90604, 90651-03, 90651A-17 (Novamatrix based)	FQ-92661320	10 / 3	7 pin male Hypertronics		5
	90419-01, 90501, 90651A-07, 90351-D/Gm; 90465; 90466	FQ-92660320	10 / 3	12 pin male Hypertronics		43
		FQ-92662320	10 / 3	7 pin male Hypertronics (must use adapter)		5
	90489, 90467, 90651A-08, 90466, 90465, 90351-D/G, 90309, 90367 (Nellcor based)	FQ-92663320	10 / 3	12 pin male Hypertronics		33


...Chart continued on following page.

SpO₂ Reusable Fingerprobes

Quick Reference Guide and Connector Chart

...Chart continued from previous page.

Original Manufacturer	Monitor Model Number	Part Number	Length ft / m	Probe Counterpart Description	Picture
Spacelabs	90351-0/6, 90465, 90466, 90467, 90489, 90651A, IM77; PC Express (Nellcor based)	FQ-92664120	3 / 0.9	DB9 male (must use adapter)	3
	PC Scout, 90303B (Nellcor based)	FQ-92665320	10 / 3	12 pin male Hypertronics	28
	90651-08, 90496 Ultraview (straight through)	FQ-92666320	10 / 3	10 pin female Spacelabs	29
	90651-08, 90496 Ultraview (Nellcor based)	FQ-92667120	3 / 0.9	DB9 male (must use adapter)	3
	90651-08, 90496 Ultraview (Novamatrix based)	FQ-92668320	10 / 3	7 pin male Hypertronics (must use adapter)	5
	90651-08, 90496 Ultraview (Novamatrix based)	FQ-92669320	10 / 3	10 pin Spacelabs female	29
S&W	All	FQ-92671320	10 / 3	10 pin female Nicolay	8
Takaoka	All (BCI based)	FQ-92681320	10 / 3	7 pin male Redel	1
		FQ-92682320	10 / 3	5 pin female Switchcraft	20
	All Monitors (BCI based)	FQ-92683120	3 / 0.9	DB9 male (must use adapter)	3
VitaLog	All (Criticare based)	FQ-92691320	10 / 3	5 pin male Lemo	6
	Vital 03	FQ-92692220	5 / 1,6	10 pin RJ45	
Weinmann	Oxycount mini (BCI based)	FQ-92701120	3 / 0.9	DB9 male	32
Welch Allyn	Propaq 206-EL, 204-EL, 202-EL, 106-EL (Nellcor based)	FQ-92711120	3 / 0.9	DB9 male	3
	62000, 652SXX, 621SO/P, 622SO/P, 623SP (Nonin based)	FQ-92713120	3 / 0.9	DB9 male	32
	62000, 652SXX, 621SO/P, 622SO/P, 623SP (Nonin based)	FQ-92714320	10 / 3	DB9 male	32
Zoll	Basic M & AED series	FQ-92721320	10 / 3	14 pin 3M plug	4



*Types

FC = Finger Clip
 FP = Finger Pediatric
 FS = Finger Silicone (Silc touch)
 PS = Pediatric Silicone (Silc touch Pedi)
 FB = Finger Silicone Boot
 FL = Finger Large Silicone Boot