

# Repair Services/Blenders



## BLENDERS

REPAIR AND OVERHAUL SERVICE PROVIDED ON:

BIRD  
IMPACT  
SECHRIST  
SIEMENS  
DATEX-OHMEDA  
PURITAN-BENNETT  
BIOMED DEVICES  
TRI-ANIM

90-DAY  
LIMITED  
WARRANTY



#09-0060  
Replacement 15 foot Oxygen  
hose, female D.I.S.S. x male  
D.I.S.S.



#09-2899  
Replacement 15 foot Medical air  
hose, female D.I.S.S. x female  
D.I.S.S.



#69-500-01  
Compressed medical air water trap  
with drain tube. Helps to remove  
moisture from medical air supply.  
Automatic operation.

# Repair Services/Regulators/Flowmeters

(Call Mercury's Biomedical  
Department for pricing)



(Call our Biomedical  
Department for pricing)

## FLOWMETERS:

AMVEX  
BIRD  
CHEMETRON  
GENSTAR  
MADA  
OHIO  
OXEQUIP  
PRECISION MEDICAL  
SP MEDICAL  
SUNMED  
TIMETER  
VACUMED  
WESTERN MEDICA



(Call our Biomedical  
Department for pricing)

## REGULATORS:

AMVEX  
AMERICAN HEALTHCARE  
B & F MEDICAL  
BIRD  
CHEMETRON  
CONCO  
ELDER  
FLOTEC  
GENSTAR  
HAMMON  
HARRIS  
HUDSON (TFX MEDICAL)  
LIFE SUPPORT PRODUCTS  
MADA  
OHIO  
OXEQUIP  
PRECISION  
PURITAN-BENNETT  
PM MEDICAL  
ROBERTSHAW  
SP MEDICAL  
TIMETER  
VICTOR  
WESTERN MEDICA



**MARK 8  
MARK 14  
MARK 25  
LARGE DIAL  
HALOSCALE COMPACT  
HALOSCALE 'INFANTA'**



**MARK 12  
MARK 20  
MARK 30  
PHYSIOLOGICAL  
HALOSCALE STANDARD  
HALOSCALE 'MIX'**

**10-DAY AVERAGE  
TURNAROUND TIME.**

**SERVICE INCLUDES:** Complete disassembly and inspection of components, cleaning, required lubrication, mercury replacement, alignment of jewels, and calibration.

## REPLACEMENT ACCESSORIES



#19-00-304  
Adapter  
23mm female  
x catheter



#19-00-410  
Adapter  
23mm female x  
15mm female



#AK-19-00-307  
Adapter  
22mm female x  
22mm female



#AK-106-7  
22mm female x  
catheter



#19-00308  
Adapter  
22mm female x  
15mm female



#AK-106-9  
22mm female x  
22mm male/15mm  
female



#19-00-409  
Adapter  
23mm male x  
23mm male



#19-00-300  
Adapter  
23mm female x  
22mm male



#19-00-481  
Adapter  
22mm female x 5/8"  
corrugated tubing

#AK-114008  
22mm male x  
15mm female

#AK-114007  
22mm female x  
15mm male

#AK-114009  
22mm male  
15mm male

## REPLACEMENT ACCESSORIES (not shown)

Case w/foam insert      #AK-108517  
Safety Harness          #AK-100S28

**NOTE:** Current model Respirometers have 22mm male inlets and outlets. However, some older models have 23mm male inlets and outlets. Older models include Mark 8, Mark 12, and physiological models. These units have a reflective chrome finished housing.

## CONTINUOUS TYPE

Oxequip  
Chemetron (NCG)  
Ohmeda  
Timeter  
Boehringer  
Precision Medical  
Western  
Amvex  
Genstar  
SunMed



## INTERMITTENT TYPE

Chemetron (NCG)  
Oxequip  
Ohmeda  
Precision Medical  
Amvex  
Genstar

*Service on above listed companies includes:  
Complete disassembly and inspection of  
components, cleaning, replacement of worn or  
defective parts and calibration.*