

Emergency Products



FLO₂™ ADJUSTABLE EMERGENCY ROOM OXYGEN SYSTEM

The FLO₂™ **BN-P60000** - 10/box is a high-concentration, high-flow oxygen diluter kit. The kit consisting of a FLO₂ Oxygen Diluter Head, a six foot respiratory hose and an oxygen mask. The FLO₂ Oxygen system allows for fast delivery of controlled oxygen concentrations from 32% to 95% with corresponding flows from 10 LPM to 41 LPM. The FLO₂ Diluter Head is also available without the hose and mask (**BN-P61000**) - 20/box for routine Respiratory Therapy use.


Oxy-PEEP™ ADJUSTABLE HIGH FLOW OXYGEN DILUTER WITH PEEP

The Oxy-PEEP **BN-P63000** is a high-concentration, high-flow oxygen diluter with adjustable PEEP capabilities. The Oxy-PEEP concept was developed by RCP's searching for a disposable, non-capital-equipment means of delivering EPAP/PEEP therapy in the emergency room and out in the field. The Oxy-PEEP system incorporates a High-Flow Oxygen Diluter Head adjustable from 32% to 95% FIO₂, 800 ml reservoir bag, a one way safety inhalation valve and an adult sealed PEEP mask with head strap adjustable from 5-20 cm of PEEP.



Asherman Chest Seal™

The ACS™ is a safer approach to treating open pneumothorax and preventing tension pneumothorax in chest injuries from gunshot, stab wounds or other penetrating chest trauma. The unique one-way valve is designed to let air and blood escape while preventing re-entry of either. Proven to be a more effective treatment than standard petroleum gauze, the ACS™ is standard issue for the U.S. Army, the British Army and the U.S. Navy.



#12-849100

10/Box