

Oxygen Hoods

SuperDome™

- Unique “one-size-fits-all” design
- Multiple oxygen ports for convenience in supplying gas to patient
- Disposable hood with reusable weights



A Solution for Oxygen Therapy in NICU

The SuperDome™ is designed to simplify oxygen therapy needs for infants and neonates. The perforated flaps allow our “one-size-fits-all” design to work in any application. The SuperDome™ is made from high-quality medical grade plastic and is held in place by two weights. We have also designed inlet holes on both sides of the dome.

These special holes will allow you to place an oxygen sensor near the face of the patient for accurate fiO2



monitoring. The SuperDome™ is also equipped with two oxygen ports located on the front of the dome providing convenient options for connecting the deflectors to the gas source.

Specifications

Dimensions (W x H x D):	11" x 6.5" x 10.25"
Weight with Weights:	0.25 lbs.
Weight without Weights:	0.18 lbs.

NICU Disposable Domes

PART NO.	MODEL	DESCRIPTION
#HO-R300P06	SuperDome	Disposable Oxygen Dome 10/pk (no weights)
#HO-R300P06-001	SuperDome Starter Kit	Disposable Oxygen Dome 10/pk (w/2 sets of weights)
#HO-R300P07	SuperDome Sample	Superdome Single pack (w/1 set of weights)

Accessories

PART NO.	MODEL	DESCRIPTION
#HO-R300P09	SuperDome Weights	5 Sets of Weights for SuperDome (10/PK)
#HO-R300P10	SuperDome Flap	Perforated Flap for SuperDome (10/Pk)
#HO-R300P11	SuperDome Deflectors	22mm Deflectors for SuperDome (10/Pk)
#HO-R300P12	SuperDome Heavy Weights	2 sets of Weights (4/PK)