


Protectors

Syphygmo-Saver




The Syphygmo-Saver was designed to prevent damage to the gauge used with blood pressure kits.
 #HH-05-250-010 Blue gauge guard
 #HH-05-250-030 Gray gauge guard
 #HH-05-250-050 Orange gauge guard

#BJ-3300 Hudson Analyzer Protector



The 3300 Hudson Analyzer Protector was designed to protect the #5577 and #6477 Hudson oxygen analyzers.

#BJ-3400 Mini-OX Analyzer Protector




The 3400 Mini-OX Analyzer protector was designed to protect the Catalyst Mini-OX I, Mini-OX II and M25RAD oxygen analyzers.




#BJ-3401
 Designed for:
 Mini-OX III Analyzer and Mini-OX IV & V Pulse Oximeters



#BJ-3402
 Designed for: Teledyne TED 60 Analyzers



#BJ-3407
 Designed for: All Fisher Paykel Nerve Stimulators



#BJ-3404
 Designed for: Hamilton Medical Max Ventilator



#BJ-3405
 Designed for: Siemens MicroO₂ and Nellcor N-30 Pulse Oximeters



#BJ-3406
 Designed for: Novametrix 510 Pulse Oximeters



#BJ-3403
 Designed for: Nonin 8500 and Healthdyne 920 Pulse Oximeters



#BJ-3408
 Designed for: Data General 5221 handheld point of service computer



#BJ-3409
 Designed for: BCI 3301 & 3302, Medical Systems Corp. 5001, Armstrong, Allied and other Pulse Oximeters.

Protectors

For Regulators



May be attached to any single-handled "D" or "E" type cart with a handle diameter of 3/4" or 1" in seconds.



The lightweight all steel gauge guards are designed to prevent the regulator/gauge from directly striking other objects. Simply connect the guard to any 3/4" or 1" tubular handle. Available in painted or chrome-plated finish.

#66-611GC-34
Green Gauge Guard for 3/4" tube handle


#66-611GC-1
Green Gauge Guard for 1" tube handle

#66-6C11GC-34
Chrome Gauge Guard for 3/4" tube handle

#66-6C11GC-1
Chrome Gauge Guard for 1" tube handle

ERIE OFFERS A BETTER, MORE ECONOMICAL WAY TO PROTECT REGULATORS


Combining the toughness of nylon with the strength and rigidity of fiberglass, ERIE'S new Gauge Guard protects delicate regulators from damage due to bumping, dropping and jarring that comes with routine use. Durable as steel, this new nylon unit flexes under pressure and reduces possible fracture of cylinder valve post. Installation is simple on most currently available regulators.



- Prevents damage and costly repairs to gauges and flowmeters.
- Attaches easily to virtually any regulator.
- Attractive green finish is easily kept clean.
- Compatible with both 1-1/2" and 2" gauge type & flowmeter regulators.

#IC-407375
Guard for regulator with 1-1/2" Gauge

#IC-407376
Guard for regulator with 2" Gauge



#AW-G901
Regard™ 1 Regulator Guard

At Timeter, we make the Regard 1 strong. For the construction of the Regard 1, we chose a heavy duty solid steel wire that is strong enough to protect the regulator by absorbing most of the shock.

The Regard 1 not only works well but looks good too. For enhanced appearance and protection from rust and corrosion, the Regard 1 is chrome-plated.

One of the best features about the Regard 1 is that it requires no tools for attaching it to the unit.

Protectors

The Gauge-Savers were designed to help eliminate unnecessary repairs due to the abuse regulators are subject to.

These 1/4" shock absorbing rubber protectors can be put on in a few seconds and DO NOT have to be removed from the regulator during cylinder changes.

Protectors are designed to fit most regulators, including:

HUDSON
VICTOR
VERIFLO
CHEMETRON
PURITAN-BENNETT

WESTERN MEDICA
MADA
TIMETER
PRECISION MEDICAL
CONCOA & MORE



Prevent this!

IMPORTANT FACT

Regulator Protectors are solely designed to minimize damage to a regulator if dropped or knocked over.

Sudden jarring or heavy impacts may cause a regulator to lose calibration or cause internal damage.

1/4" shock absorbing rubber

BJ-2901 Gauge Saver



Fits 1 1/2" BACK port gauges with twist-on lens.

BJ-2911 Gauge Saver



Fits 1 1/2" BACK port gauges with snap-on lens.

BJ-2904 Gauge Saver



Fits 1 1/2" BOTTOM port gauges with twist-on lens.

BJ-2914 Gauge Saver



Fits 1 1/2" BOTTOM port gauges with snap-on lens.

BJ-2909 Gauge Saver



Fits all 2" BOTTOM port gauges.

BJ-2902 Gauge Saver




Fits both 2" side by side gauges

BJ-2903 Gauge Saver



Fits both 1 1/2" side by side gauges with twist on lens

BJ-2903-W Gauge Saver



Fits both 1 1/2" side by side gauges with snap on lens

BJ-2905 Knob Protector

Fits 1 3/4" control knobs.



Western Enterprise
Veriflo Hudson
"newer" 2600
Series and more

BJ-2915 Knob Protector

Fits 2" control knobs.



Victor "newer"
VMG Series
Chemetron
precision
medical and more.

Protectors


BJ-2800 The Gage-Gard was designed to help eliminate accidental damage to manometers. The Gage-Gard is supplied with a Velcro type of fastener for temporary or permanent mounting to any smooth Surface.

Model #BJ-2800

Most 2-1/2" BOTTOM port pressure manometers, OEM, DHD, Portex, Instrumentation Industries, Anesthesia Associates and many others.



#BJ-2800



#BJ-2915
Fits 2" control knobs.
Victor "newer" VMG Series,
Chemetron and more.



#BJ-2801 Gage-Gard
Most 2-1/2" Back port pressure manometers, Bird, Monaghan, Puritan-Bennett Instrumentation Industries and many others.

#34-G-15A-PCE
Gauge Saver designed for 1-1/2" port Slim-line Gauges.

#BJ-2801



#BJ-3000

The Bo-Re-Guard was designed for all Boehringer Spirometers. A lanyard is also included.

#BJ-2909
Fits all 2" BOTTOM port gauges.



#BJ-2909

#BJ-3000 Bo-Re-Guard

The BJ-2703 carrying case was specifically designed Extra Large for the purpose of holding respirometers with Spiro-Guards on them. The plastic case does NOT have a zipper and is lightweight and durable.

(Respirometer, Spiro-Guard and adapters not included.)

9" x 5" with carrying case.



#BJ-2703

Protectors

Respirometers



#BJ-2700 Spiro-Guard

The Spiro-Guard was specifically designed to help eliminate accidental damage to mechanical respirometers that can occur in the hurried environment of intensive care. This 1/4" shock absorbing rubber is very effective and will help to keep your respirometers "In Service" rather than "Out of Service."



#BJ-2701 Spiro Guard

MODEL#

DESIGNED FOR:

#BJ-2700

Mark 8, Mark 14, Mark 25,
Physiological, Haloscale Compact

#BJ-2701

Fraser-Harlake RM 121, 311
(Ohmeda)

#BJ-2702

Mark 12, Mark 20,
Haloscale Standard



#BJ-2702 Spiro Guard

Each Spiro-Guard is shipped with rubber O-rings and a lanyard.



#BJ-270025 Universal Adapter

Combination of BJ-270030 and BJ-270031



#BJ-270030 Polyurethane Adapter

22 - 23mm x 22 - 23mm Accommodates either 22mm or 23mm outlets



#BJ-270031 Silicone Adapter

24mm x 15mm